

**Reporting Title:** Huntington Disease Analysis

**Performing Location:** Rochester

**Shipping Instructions:**

Specimen preferred to arrive within 96 hours of collection.

**Necessary Information:**

[Molecular Genetics: Neurology Patient Information](#) or a recent clinical note is required. Testing cannot proceed without this information.

**Specimen Requirements:**

**Patient Preparation:** A previous bone marrow transplant from an allogenic donor will interfere with testing. Call 800-533-1710 for instructions for testing patients who have received a bone marrow transplant.

**Specimen Type:** Whole blood

**Container/Tube:**

**Preferred:** Lavender top (EDTA) or yellow top (ACD)

**Acceptable:** Any anticoagulant

**Specimen Volume:** 3 mL

**Collection Instructions:**

1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. **Do not aliquot.**

**Forms:**

1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available:

- [Informed Consent for Genetic Testing](#) (T576)
- [Informed Consent for Genetic Testing-Spanish](#) (T826)

2. If not ordering electronically, complete, print, and send a [Neurology Specialty Testing Client Test Request](#) (T732) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Frozen		
	Refrigerated		

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
52885	Result Summary	Alphanumeric		50397-9
52886	Result	Alphanumeric		53782-9
52887	Interpretation	Alphanumeric		69047-9
52888	Reason for Referral	Alphanumeric		42349-1
52889	Specimen	Alphanumeric		31208-2
52890	Source	Alphanumeric		31208-2
52891	Released By	Alphanumeric		18771-6

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

81271-HTT (huntingtin) (eg, Huntington disease) gene analysis; evaluation to detect abnormal (eg, expanded) alleles

**Reference Values:**

Normal alleles: <27 CAG repeats

Intermediate alleles: 27-35 CAG repeats

Reduced penetrance: 36-39 CAG repeats

Full penetrance: >39 CAG repeats

An interpretive report will be provided.