

Test Definition: BAIPD

Bile Acids for Peroxisomal Disorders, Serum

Reporting Title: Bile Acids for Peroxisomal D/O, S

Performing Location: Rochester

Ordering Guidance:

For assessment of general liver dysfunction in adults or diagnosis or monitoring of intrahepatic cholestasis of pregnancy, order BAFS / Bile Acids, Fractionated and Total, Serum.

Specimen Requirements:

Patient Preparation: Patient must be fasting for 12 to 14 hours.

Collection Container/Tube:

Preferred: Serum gel **Acceptable:** Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL

Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

Forms:

If not ordering electronically, complete, print, and send a Biochemical Genetics Test Request (T798) with the specimen.

Specimen Type	Temperature	Time	Special Container	
Serum	Refrigerated (preferred)	90 days		
	Ambient	90 days		
	Frozen	90 days		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
41446	Dihydroxycholestanoic Acid	Numeric	nmol/mL	53479-2
41447	Trihydroxycholestanoic Acid	Numeric	nmol/mL	38188-9
41448	Total Cholic Acid	Alphanumeric	nmol/mL	30518-5
41449	Total Chenodeoxycholic Acid	Alphanumeric	nmol/mL	30519-3
41450	Total Ursodeoxycholic Acid	Alphanumeric	nmol/mL	55159-8
41451	Total Bile Acids	Alphanumeric	nmol/mL	14628-2

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82542

Reference Values:

Dihydroxycholestanoic acid: < or =0.10 nmol/mL



Total bile acids: < or =19.00 nmol/mL

Test Definition: BAIPD

Bile Acids for Peroxisomal Disorders, Serum

Trihydroxycholestanoic acid: < or =1.30 nmol/mL
Total cholic acid: < or =5.00 nmol/mL
Total chenodeoxycholic acid: < or =6.00 nmol/mL
Total ursodeoxycholic acid: < or =2.00 nmol/mL