

Reporting Title: Pyridoxic Acid (PA), P
Performing Location: Rochester

Shipping Instructions:
Ship specimen in amber vial to protect from light.

Specimen Requirements:
Only orderable as part of a profile. For more information see B6PRO / Vitamin B6 Profile (Pyridoxal 5-Phosphate and Pyridoxic Acid), Plasma.

Patient Preparation:
1. Fasting-overnight (12-14 hours) (infants-collect specimen prior to next feeding). Water can be taken as needed.
2. For 24 hours prior to specimen collection, patient must not take multivitamins or vitamin supplements.

Supplies: Amber Frosted Tube, 5 mL (T915)
Collection Container/Tube: Green top (sodium or lithium heparin) or plasma gel separator (PST)
Submission Container/Tube: Amber vial
Specimen Volume: 1 mL

Collection Instructions:
1. Centrifuge at 4 degrees C within 2 hours of collection.
2. Aliquot all plasma into amber vial and freeze immediately.

Specimen Type	Temperature	Time	Special Container
Plasma Heparin	Frozen	29 days	LIGHT PROTECTED

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
61065	Pyridoxic Acid (PA), P	Numeric	mcg/L	1688-1

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
82542

Reference Values:
Only orderable as part of a profile. For more information see B6PRO / Vitamin B6 Profile (Pyridoxal 5-Phosphate and Pyridoxic Acid), Plasma.

3-30 mcg/L