

11-Deoxycortisol, Serum

# **Reporting Title:** 11-Deoxycortisol, S **Performing Location:** Rochester

## **Necessary Information:**

Indicate if specimen was collection before or after metyrapone administration.

### **Specimen Requirements:**

Supplies: Sarstedt Aliquot Tube, 5 mL (T914) Collection Container/Tube: Preferred: Red top Acceptable: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL

#### **Collection Instructions:**

- 1. Morning (8 a.m.) specimen is preferred.
- 2. Centrifuge and aliquot serum into a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

## Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
46923	11-Deoxycortisol, S	Numeric	ng/dL	1657-6

LOINC® and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

**CPT Code Information:** 82634

Reference Values: < or =18 years: <344 ng/dL >18 years: 10-79 ng/dL

For International System of Units (SI) conversion for Reference Values, see www.mayocliniclabs.com/order-tests/si-unit-conversion.html.