

Hepatitis B Virus Genotyping

## **Reporting Title:** Hepatitis B Virus Genotype **Performing Location:** ARUP Laboratories

## **Specimen Requirements:**

Draw blood in lavender (EDTA) tube(s). Spin down and send 2 mL plasma frozen in a plastic vial.

### Required: 1. Viral Load

2. Viral Load Date

**Note**: Red-top serum and serum gel tube(s) are acceptable.

Note: This test may be unsuccessful if the HBV Viral load is less than log 3.0 or 1,000 IU/mL of plasma.

Specimen Type	Temperature	Time	Special Container
Plasma EDTA	Frozen (preferred)	42 days	
	Refrigerated	7 days	
	Ambient	72 hours	

## Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
FHBG	Z3193	HBV Viral Load	Plain Text	Yes
FHBG	Z3194	HBV Viral Load Date	Plain Text	Yes

# Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3195	Hepatitis B Genotype	Alphanumeric		32366-7
Z3196	HBV Surface Antigen Mutations	Alphanumeric		32366-7
Z3197	HBV RT Polymerase Mutations	Alphanumeric		Unable to Verify

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

## **CPT Code Information:**

87912

#### **Reference Values:**

Interpretive Information: Hepatitis B Virus Genotype

HBV genotype and resistance interpretation is provided by SeqHepB software from Evivar Medical. The following



# **Test Definition: FHBG**

Hepatitis B Virus Genotyping

mutations are reported: reverse transcriptase L80I/V, II69T, V173L, L180M, A181S/T/V, T184A/C/F/I/G/S/M/L, S202C/G/I, M204I/V, N236T, M250I/L/V; surface antigen P120T, D144A, G145R.

Both the HBV RT polymerase and the HBsAg encoding regions are sequenced. Resistance and surface antigen mutations are reported. In addition, the major HBV genotypes are identified. Mutations in viral sub-populations below 20% of total may not be detected.

This test is performed pursuant to a license agreement with Roche Molecular Systems, Inc.