

Test Definition: FZOLP

Zolpidem (Ambien), serum or plasma

Reporting Title: Zolpidem (Ambien), Serum **Performing Location:** Medtox Laboratories, Inc.

Specimen Requirements:

Submit only 1 of the following specimens:

Plasma

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 2 mL sodium heparin plasma refrigerated in a plastic vial.

Serum

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Frozen	180 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3322	Zolpidem, SE	Alphanumeric		6953-4

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80368

Reference Values:

Units: ng/mL

Expected hypnotic zolpidem concentrations in patients taking recommended daily dosages: up to 250 ng/mL.

Toxic range has not been established.