

Test Definition: VWACT

von Willebrand Factor Activity, Plasma

Reporting Title: von Willebrand Factor Activity, P

Performing Location: Rochester

Ordering Guidance:

This activity assay is most effective when it is combined with measurement of von Willebrand factor antigen and factor VIII coagulant activity, preferably as a panel of tests with reflexive testing and interpretive reporting. See AVWPR / von Willebrand Disease Profile, Plasma.

Additional Testing Requirements:

Tests for F8A / Coagulation Factor VIII Activity Assay, Plasma and VWAG / von Willebrand Factor Antigen, Plasma are recommended in conjunction with this test (von Willebrand activity).

Specimen Requirements:

Specimen Type: Platelet-poor plasma

Collection Container/Tube: Light-blue top (3.2% sodium citrate) **Submission Container/Tube:** Plastic vial (polypropylene preferred)

Specimen Volume: 2 mL in 2 vials each containing 1 mL

Collection Instructions:

- 1. Specimen must be collected prior to factor replacement therapy.
- 2. For complete instructions, see Coagulation Guidelines for Specimen Handling and Processing.
- 3. Centrifuge, transfer all plasma into a polypropylene vial, and centrifuge plasma again.
- 4. Aliquot plasma (1-2 mL per aliquot) into 2 separate polypropylene vials leaving 0.25 mL in the bottom of centrifuged vial.
- 5. Freeze plasma immediately (no longer than 4 hours after collection) at -20 degrees C or, ideally, < or =-40 degrees C.
- 6. Send specimens in the same shipping container.

Additional Information:

- 1. Double-centrifuged specimen is critical for accurate results as platelet contamination may cause spurious results.
- 2. Each coagulation assay requested should have its own vial.

Forms:

If not ordering electronically, complete, print, and send a Coagulation Test Request (T753) with the specimen.

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

Result Codes:

1	Result ID	Reporting Name	Туре	Unit	LOINC®
١	VWACT	von Willebrand Factor Activity, P	Numeric	%	68324-3

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No



Test Definition: VWACT

von Willebrand Factor Activity, Plasma

CPT Code Information:

85397

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
RIST	Ristocetin Cofactor, P	1	85245	No	No

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
RIST	9046	Ristocetin Cofactor, P	Numeric	%	6014-5

Reference Values:

55-200%

Normal, full-term newborn infants may have mildly increased levels which reach adult levels by 90 days postnatal. Healthy, premature infants (30-36 weeks gestation) may have increased levels that reach adult levels by 180 days.

Note: Individuals of blood group "O" may have lower plasma von Willebrand factor (VWF) activity than those of other ABO blood groups, such that apparently normal individuals of blood group "O" may have plasma VWF activity as low as 40% to 50%, whereas the lower limit of the reference range for individuals of other blood groups may be 60% to 70%.