

# **Test Definition: DRVI2**

Dilute Russell's Viper Venom Time (DRVVT)

Mix Ratio, Plasma

**Reporting Title:** DRVVT Mix Ratio **Performing Location:** Rochester

# **Specimen Requirements:**

Only orderable as part of a reflex. For more information see DRVI1 / Dilute Russell's Viper Venom Time (DRVVT), with Reflex, Plasma.

Specimen Type: Platelet-poor plasma

**Collection Container/Tube:** Light-blue top (3.2% sodium citrate)

Submission Container/Tube: Plastic vial

**Specimen Volume:** 1 mL **Collection Instructions:** 

- 1. For complete instructions, see Coagulation Guidelines for Specimen Handling and Processing.
- 2. Centrifuge, remove plasma, and centrifuge plasma again.
- 3. Aliquot into a separate plastic vial, leaving 0.25 mL in the bottom of the centrifuged vial.
- 4. Freeze plasma immediately (no longer than 4 hours after collection) at -20 degrees C or, ideally, -40 degrees C or below.

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
RMRI2	DRVVT Mix Ratio	Numeric	ratio	75512-4

LOINC® and CPT codes are provided by the performing laboratory.

# Supplemental Report:

No

## **CPT Code Information:**

85613

#### **Reference Values:**

Only orderable as part of a reflex. For more information see DRVI1 / Dilute Russell's Viper Venom Time (DRVVT), with Reflex, Plasma.

#### <1.20

Normal ranges for children: not clearly established, but similar to normal ranges for adults, except for newborn infants whose results may not reach adult values until 3 to 6 months of age.