

**Reporting Title:** PLA2R, Immunofluorescence, S  
**Performing Location:** Rochester

**Ordering Guidance:**  
This test can be used in a complementary fashion with a quantitative enzyme-linked immunosorbent assay to confirm a positive result, especially if it is a low-level titer. Since the immunofluorescence is more sensitive, it can also be used to detect or follow patients with low level antibodies not detected by ELISA.

**Specimen Requirements:**  
**Supplies:** Sarstedt Aliquot Tube 5 mL (T914)  
**Collection Container/Tube:**  
**Preferred:** Serum gel  
**Acceptable:** Red top  
**Submission Container/Tube:** Plastic vial  
**Specimen Volume:** 1 mL  
**Collection Instructions:** Centrifuge and aliquot serum into a plastic vial within 2 hours of collection.

**Forms:**  
[If not ordering electronically, complete, print, and send a Renal Diagnostics Test Request](#) (T830) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	8 hours	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
PLA2I	PLA2R, Immunofluorescence, S	Alphanumeric		82991-1

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**  
No

**CPT Code Information:**  
86255

**Reference Values:**  
Negative