

# **Test Definition: 11INE**

Factor XI Inhibitor Evaluation, Plasma

Reporting Title: Factor XI Inhib Profile, P

Performing Location: Rochester

### **Ordering Guidance:**

This test is for factor XI inhibitors only. If the presence or type of inhibitor is unknown, first order APROL / Prolonged Clot Time Profile, Plasma, except for patients with known hemophilia A or B. When screening studies are needed for patients with known hemophilia A or B, order 8INHE / Factor IX Inhibitor Evaluation, Plasma; or 9 INHE / Factor IX Inhibitor Evaluation, Plasma; respectively.

### **Shipping Instructions:**

Send all vials in the same shipping container.

### **Necessary Information:**

If priority specimen, mark request form, give reason, and request a call-back.

### **Specimen Requirements:**

Specimen Type: Platelet-poor plasma

### **Patient Preparation:**

1. Patient should not be receiving Coumadin (warfarin), heparin, direct thrombin inhibitors (argatroban, dabigatran), or direct factor Xa inhibitors (apixaban, rivaroxaban, and edoxaban).

2. Fasting preferred.

**Collection Container/Tube:** Light-blue top (3.2% sodium citrate)

Submission Container/Tube: 3 Plastic vials

Specimen Volume: 3 mL in 3 plastic vials, each containing 1 mL

#### **Collection Instructions:**

- 1. Specimen must be collected prior to factor replacement therapy.
- 2. For complete instructions, see Coagulation Guidelines for Specimen Handling and Processing.
- 3. Centrifuge, transfer all plasma into a plastic vial, and centrifuge plasma again.
- 4. Aliquot plasma (1-2 mL per aliquot) into 3 separate plastic vials, leaving 0.25 mL in the bottom of centrifuged vial.
- 5. Freeze plasma immediately at -20 degrees C or, ideally, at -40 degrees C or below.

### **Additional Information:**

- 1. Double-centrifuged specimen is critical for accurate results as platelet contamination may cause spurious results.
- 2. Each coagulation assay requested should have its own vial.

#### Forms:

If not ordering electronically, complete, print, and send a Coagulation Test Request (T753) with the specimen.

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

# **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
F_11	Coag Factor XI Assay, P	Numeric	%	3226-8
11INT	FXI Inhib Profile Tech Interp	Alphanumeric		69049-5



# **Test Definition: 11INE**

Factor XI Inhibitor Evaluation, Plasma

LOINC® and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

## **Components:**

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
11INT	FXI Inhib Profile Tech Interp	1	85390	Yes	No
F_11	Coag Factor XI Assay, P	1	85270	Yes	Yes

# **CPT Code Information:**

85390

85270

85335 (if appropriate)

85335 (if appropriate)

85390 (if appropriate)

### **Reflex Tests:**

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
11AIH	FXI Inhib Profile Prof Interp	1	85390	No	No
11_IS	Factor XI Inhib Scrn	1	85335	No	No
GBETH	General Factor Bethesda Units, P	1	85335	No	No

## **Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
11_IS	7804	Factor XI Inhib Scrn	Alphanumeric		80603-4
11AIH	607492	Reviewed by	Alphanumeric		18771-6
11AIH	607448	FXI Inhib Profile Prof Interp	Alphanumeric		69049-5
GBETH	607434	General Factor Bethesda Units, P	Numeric	BU	13591-3

# **Reference Values:**

**FACTOR XI ACTIVITY ASSAY** 

Adults: 55-150%

Normal, full-term newborn infants or healthy premature infants may have decreased levels (> or =10%) that may not reach adult levels for 180 days or more postnatal.\*

### FACTOR XI INHIBITOR SCREEN:

Negative

<sup>\*</sup>See Pediatric Hemostasis References section in Coagulation Guidelines for Specimen Handling and Processing.



# **Test Definition: 11INE**

Factor XI Inhibitor Evaluation, Plasma

GENERAL FACTOR BETHESDA UNITS: < or =0.5 Bethesda Units