

**Reporting Title:** Arsenic/Creat w/Reflex, Random,U  
**Performing Location:** Rochester

**Specimen Requirements:**

**Patient Preparation:**

- Patient should not eat seafood for a 48-hour period prior to start of collection.
- High concentrations of gadolinium and iodine are known to interfere with most metal tests. If either gadolinium- or iodine-containing contrast media has been administered, a specimen should not be collected for 96 hours.

**Supplies:** Urine Tubes, 10 mL (T068)

**Collection Container/Tube:** Clean, plastic urine container with no metal cap or glued insert

**Submission Container/Tube:** Plastic, 10-mL urine tube or clean, plastic aliquot container with no metal cap or glued insert

**Specimen Volume:** 6 mL

**Collection Instructions:**

1. Collect urine a random urine specimen.
2. See [Metals Analysis Specimen Collection and Transport](#) for complete instructions.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	28 days	
	Frozen	28 days	
	Ambient	72 hours	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
CRETR	Creatinine, Random, U	Numeric	mg/dL	2161-8
608900	Arsenic/Creatinine Ratio, U	Numeric	mcg/g Cr	13463-5
608901	Total Arsenic Concentration	Numeric	mcg/L	5586-3

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**Components:**

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
ASCU	Arsenic/Creatinine Ratio, U			Yes	No
CRETR	Creatinine, Random, U	1	82570	Yes	No

**CPT Code Information:**

82175

82570

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
SPAS	Arsenic Speciation, Random, U	1	82175	No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
SPAS	607692	Inorganic Arsenic (Toxic)	Alphanumeric	mcg/L	12481-8
SPAS	607693	Methylated Arsenic (Toxic)	Alphanumeric	mcg/L	53779-5
SPAS	607694	Organic Arsenic (Non-Toxic)	Alphanumeric	mcg/L	53778-7
SPAS	609405	Toxic Arsenic Concentration	Alphanumeric	mcg/L	96250-6
SPAS	607691	Total Arsenic Concentration	Numeric	mcg/L	5586-3

Reference Values:

0-17 years: Not established  
> or =18 years: <24 mcg/g creatinine