
Reporting Title: CYP2C9 Genotype, V

Performing Location: Rochester

Ordering Guidance:

If patient is or will be using warfarin, the preferred test is WARSQ / Warfarin Response Genotype, Varies, which includes testing of *CYP2C9*, *VKORC1*, *CYP4A2*, and rs12777823.

Testing is available as the single gene assay (this test) or as a part of a focused pharmacogenomics panel, which includes testing for the following genes: *CYPs 1A2, 2C9, 2C19, 2D6, 3A4, 3A5, 4F2, SLCO1B1*, and *VKORC1*.

If multiple pharmacogenomic genotype testing is needed, order PGXQP / Focused Pharmacogenomics Panel, Varies.

Specimen Requirements:

Multiple genotype tests can be performed on a single specimen after a single extraction. See [Multiple Genotype Test List](#) for a list of tests that can be ordered together.

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube: Lavender top (EDTA)

Specimen Volume: 3 mL

Collection Instructions:

1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. **Do not aliquot.**

Specimen Stability Information: Ambient (preferred) 9 days/Refrigerated 30 days

Specimen Type: Saliva

Patient Preparation: Patient should not eat, drink, smoke, or chew gum 30 minutes prior to collection.

Supplies: Saliva Swab Collection Kit (T786)

Specimen Volume: 1 Swab

Collection Instructions: Collect and send specimen per kit instructions.

Specimen Stability Information: Ambient 30 days

Specimen Type: Extracted DNA

Container/Tube: 2-mL screw top tube

Specimen Volume: 100 mcL (microliters)

Collection Instructions:

1. The preferred volume is 100 mcL at a concentration of 50 ng/mcL.
2. Provide concentration of DNA and volume on tube.

Specimen Stability Information: Frozen (preferred)/Ambient/Refrigerated

Forms:

1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available:

-[Informed Consent for Genetic Testing](#) (T576)

-[Informed Consent for Genetic Testing-Spanish](#) (T826)

2. If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

- [Neurology Specialty Testing Client Test Request](#) (T732)
- [Therapeutics Test Request](#) (T831)

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
610096	CYP2C9 Genotype	Alphanumeric		46724-1
610568	CYP2C9 Activity Score	Alphanumeric		In Process
610098	Interpretation	Alphanumeric		69047-9
610099	Additional Information	Alphanumeric		48767-8
610100	Method	Alphanumeric		85069-3
610101	Disclaimer	Alphanumeric		62364-5
610102	Reviewed by	Alphanumeric		18771-6
610097	CYP2C9 Phenotype	Alphanumeric		79716-7

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

81227

Reference Values:

An interpretive report will be provided.