

Test Definition: CRT24

Creatinine, 24 Hour, Urine

Reporting Title: Creatinine, 24 HR, U

Performing Location: Rochester

Specimen Requirements:

Only orderable as part of a profile. For more information see:

- -NMH24 / N-Methylhistamine, 24 Hour, Urine
- -RB24 / Retinol-Binding Protein, 24 Hour, Urine
- -A124 / Alpha-1-Microglobulin, 24 Hour, Urine

Urine Preservative Collection Options:

Note: The addition of preservative or application of temperature controls **must occur within 4 hours of completion** of the collection.

Ambient	Yes
Refrigerate	Preferred
Frozen	Yes
50% Acetic Acid	Yes
Boric Acid	Yes
Diazolidinyl Urea	Yes
6M Hydrochloric Acid	Yes
6M Nitric Acid	Yes
Sodium Carbonate	Yes
Thymol	No
Toluene	No

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	30 days	
	Frozen	30 days	
	Ambient	14 days	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
CRT24	TM10	Collection Duration	Plain Text	Yes
CRT24	VL8	Urine Volume	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
TM10	Collection Duration	Numeric	h	13362-9
VL8	Urine Volume	Numeric	mL	3167-4
CR_A	Creatinine, 24 HR, U	Numeric	mg/24 h	2162-6
CR_24	Creatinine Concentration, 24 HR, U	Numeric	mg/dL	20624-3



Test Definition: CRT24

Creatinine, 24 Hour, Urine

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

Reference Values:

Only orderable as part of a profile. For more information see:

- -NMH24 / N-Methylhistamine, 24 Hour, Urine
- -RB24 / Retinol-Binding Protein, 24 Hour, Urine
- -A124 / Alpha-1-Microglobulin, 24 Hour, Urine

Normal values mg per 24 hours: Males: 930-2955 mg/24 hours Females: 603-1783 mg/24 hours

Reference values have not been established for patients who are less than 18 years of age.

For SI unit Reference Values, see https://www.mayocliniclabs.com/order-tests/si-unit-conversion.html