

# **Reporting Title:** Magnesium, 24 HR, U **Performing Location:** Rochester

### **Necessary Information:**

Specimen volume in milliliters and duration are required.

## **Specimen Requirements:**

Supplies: Sarstedt 5 mL Aliquot Tube (T914)

**Collection Container/Tube:** 24-hour graduated urine container with no metal cap or glued insert **Submission Container/Tube:** Plastic, 5 mL tube or a clean, plastic urine container with no metal cap or glued insert

Specimen Volume: 4 mL

**Collection Instructions:** 

1. Collect urine for 24 hours.

2. Refrigerate specimen within 4 hours of completion of 24-hour collection.

Additional Information: See <u>Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens</u> for multiple collections.

### Forms:

If not ordering electronically, complete, print, and send a <u>Renal Diagnostics Test Request</u> (T830) with the specimen.

## **Urine Preservative Collection Options:**

**Note:** The addition of preservative or application of temperature controls **must occur within 4 hours of completion** of the collection.

Ambient	ОК
Refrigerate	Preferred
Frozen	ОК
50% Acetic Acid	ОК
Boric Acid	ОК
Diazolidinyl Urea	ОК
6M Hydrochloric Acid	ОК
6M Nitric Acid	No
Sodium Carbonate	No
Thymol	ОК
Toluene	No

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	30 days	
	Ambient	72 hours	

# Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
MAGU	TM115	Collection Duration	Plain Text	Yes
MAGU	VL111	Urine Volume	Plain Text	Yes



## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
MG24	Magnesium, 24 HR, U	Numeric	mg/24 h	24447-5
TM115	Collection Duration	Numeric	h	13362-9
VL111	Urine Volume	Numeric	mL	3167-4

LOINC® and CPT codes are provided by the performing laboratory.

### Supplemental Report:

No

### **CPT Code Information:**

83735

## **Reference Values:**

51-269 mg/24 hours

Reference values have not been established for patients <18 years and >83 years of age.