

## **Test Definition: WASEQ**

Alpha Globin Gene Sequencing, Varies

Reporting Title: Alpha Globin Gene Sequencing, B

Performing Location: Rochester

### Ordering Guidance:

For first-tier testing for alpha thalassemia detection, order THEV1 / Thalassemia and Hemoglobinopathy Evaluation, Serum and Whole Blood.

For first-tier testing for an alpha-globin variant, order HBEL1 / Hemoglobin Electrophoresis Evaluation, Blood.

If genetic testing is desired, the first-tier genetic test assesses large deletional alpha-thalassemia alterations. Order ATHAL / Alpha-Globin Gene Analysis, Varies.

### **Necessary Information:**

- 1. Patient's age is required.
- 2. Include recent transfusion information.

### **Specimen Requirements:**

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube:

**Preferred:** Lavender top (EDTA)

Acceptable: Yellow top (ACD), green top (sodium heparin)

Specimen Volume: 4 mL Collection Instructions:

1. Invert several times to mix blood.

2. Send whole blood specimen in the original tube. Do not aliquot

Specimen Stability Information: Refrigerate 30 days(preferred)/Ambient 14 days

Specimen Type: Extracted DNA from whole blood

Container/Tube: 1.5- to 2-mL tube Specimen Volume: Entire specimen

**Collection Instructions:** 

- 1. Label specimen as extracted DNA and source of specimen
- 2. Provide volume and concentration of the DNA

Specimen Stability Information: Frozen (preferred)/Refrigerate/Ambient

### Forms:

- 1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available in Special Instructions:
- -Informed Consent for Genetic Testing (T576)
- -Informed Consent for Genetic Testing-Spanish (T826)
- 2. Metabolic Hematology Patient Information (T810)
- 3. If not ordering electronically, complete, print, and send <u>Benign Hematology Test Request Form</u> (T755) with the specimen



# **Test Definition: WASEQ**

Alpha Globin Gene Sequencing, Varies

Specimen Type	Temperature	Time	Special Container	
Varies	Varies			

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
43921	Interpretation	Alphanumeric		69047-9
61362	Alpha Globin Gene Sequencing, B	Alphanumeric		87730-8

LOINC® and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

## **CPT Code Information:**

81259-HBA1/HBA2; full sequence

## **Reference Values:**

An interpretive report will be provided.