

Reporting Title: Factor H Autoantibody, S

Performing Location: Rochester

Specimen Requirements:

Collection Container/Tube:

Preferred: Red top

Acceptable: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL

Collection Instructions:

1. Immediately after specimen collection, place the tube on wet ice.
2. Centrifuge and aliquot serum into plastic vial.
3. Freeze specimen within 30 minutes.

Additional Information: If the specimen is to be shared with AHUSD / Atypical Hemolytic Uremic Syndrome Complement Panel, Serum and Plasma, only serum collected in a red-top tube is acceptable.

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------|---------|-------------------|
| Serum Red | Frozen (preferred) | 28 days | |
| | Refrigerated | 28 days | |
| | Ambient | 14 days | |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|--------------------------|---------|------|----------|
| AFH | Factor H Autoantibody, S | Numeric | U/mL | 101863-9 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83520

Reference Values:

<15.8 U/mL