

Reporting Title: Calcitonin, FNAB, Lymph Node
Performing Location: Rochester

Shipping Instructions:
Send specimen frozen.

Necessary Information:
The biopsied site of each specimen **is required** and **must be** clearly identified in the LIS and/or batch sheet.

- Specimen Requirements:**
- Patient Preparation:** For 12 hours before this procedure do not take multivitamins or dietary supplements containing biotin (vitamin B7), which is commonly found in hair, skin, and nail supplements and multivitamins.
- Container/Tube:** Plain, plastic, screw-top tube
- Specimen Volume:** 1 to 1.5 mL
- Collection Instructions:**
1. Needle wash specimens for analysis should be collected in conjunction with cytology specimens.
 2. Have saline available prior to start of procedure. Saline is the only acceptable solution for needle washings.
 3. After each fine-needle aspiration biopsy (FNAB) has been collected and the material in the needle has been expelled onto a slide for cytologic analysis, attach the used FNAB needle to an empty syringe.
 4. Withdraw between 0.10 mL and 0.25 mL of saline up through the needle until the saline starts to fill the hub of the needle or end of the syringe.
 5. Expel this fluid back through the needle into a separate tube. This is the needle washing used for analysis.
 6. Repeat steps 2 through 4 for each needle pass of the same biopsied site and empty into the same tube, accumulating a total of 0.5 mL to 1.5 mL of fluid to send to the laboratory. (If more than 1 site is biopsied, see Additional Information)
 7. Inspect specimen for visible blood or tissue contamination:
 - a. If bloody, centrifuge specimen and transfer supernatant to a new plastic aliquot tube (5-mL standard tube) to send to laboratory. The supernatant, not the cellular material, is used for analysis.
 - b. If specimen is clear, centrifugation is not necessary.
 8. Freeze within 2 to 4 hours of collection.

- Additional Information:**
1. If more than 1 site is biopsied, each washing material should be submitted on a separate tube and under a different order number.
 2. A minimum of 0.5 mL is required for testing; however, the total collection volume should not exceed 1.5 mL. Sample volumes outside these parameters may be rejected.
 3. Do not send saline control. This test has been validated to rule-out saline matrix effect.

Forms:
If not ordering electronically, complete, print, and send an [Oncology Test Request](#) (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container
Fine Needle Wash	Frozen (preferred)	7 days	
	Refrigerated	4 hours	

Ask at Order Entry (AOE) Questions:

Test Definition: CATLN

Calcitonin, Fine-Needle Aspiration Biopsy
Needle Wash, Lymph Node

Test ID	Question ID	Description	Type	Reportable
CATLN	SITEB	Site	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
CATL	Calcitonin, FNAB, Lymph Node	Numeric	pg/mL	75709-6
SITEB	Site	Alphanumeric		39111-0

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82308

Reference Values:

An interpretive report will be provided.