

## **Test Definition: CAORG**

Comprehensive Marfan, Loeys-Dietz, Ehlers-Danlos, and Aortopathy Gene Panel,

Varies

Reporting Title: Comprehensive Aortopathy Gene Panel

Performing Location: Rochester

### **Ordering Guidance:**

Customization of this panel and single gene analysis for any gene present on this panel are available. For more information see CGPH / Custom Gene Panel, Hereditary, Next-Generation Sequencing, Varies.

Targeted testing for familial variants (also called site-specific or known mutations testing) is available for the genes on this panel. See FMTT / Familial Mutation, Targeted Testing, Varies. To obtain more information about this testing option, call 800-533-1710.

### **Shipping Instructions:**

Specimen preferred to arrive within 96 hours of collection.

### **Necessary Information:**

<u>Prior Authorization</u> is available, **but not required**, for this test. If proceeding with the prior authorization process, submit the required form with the specimen.

### **Specimen Requirements:**

**Patient Preparation:** A previous bone marrow transplant from an allogenic donor will interfere with testing. Call 800-533-1710 for instructions for testing patients who have received a bone marrow transplant.

Specimen Type: Whole blood

Container/Tube:

**Preferred:** Lavender top (EDTA) or yellow top (ACD)

Acceptable: Any anticoagulant Specimen Volume: 3 mL Collection Instructions:

1. Invert several times to mix blood

2. Send whole blood specimen in original tube. **Do not aliquot. Specimen Stability Information:** Ambient (preferred)/Refrigerated

#### Forms:

- 1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available:
- -Informed Consent for Genetic Testing (T576)
- -Informed Consent for Genetic Testing (Spanish) (T826)
- 2. Connective Tissue/Cerebrovascular Disease Genetic Testing Patient Information
- 3. Comprehensive Aortopathy Gene Panel (CAORG) Prior Authorization Ordering Instructions
- 4. If not ordering electronically, complete, print, and send a <u>Cardiovascular Test Request</u> (T724) with the specimen.

Specimen Type	Temperature	Time	Special Container	
Varies	Varies			



# **Test Definition: CAORG**

Comprehensive Marfan, Loeys-Dietz, Ehlers-Danlos, and Aortopathy Gene Panel, Varies

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
617156	Test Description	Alphanumeric		62364-5
617157	Specimen	Alphanumeric		31208-2
617158	Source	Alphanumeric		31208-2
617159	Result Summary	Alphanumeric		50397-9
617160	Result	Alphanumeric		82939-0
617161	Interpretation	Alphanumeric		69047-9
617162	Additional Results	Alphanumeric		82939-0
617163	Resources	Alphanumeric		99622-3
617164	Additional Information	Alphanumeric		48767-8
617165	Method	Alphanumeric		85069-3
617166	Genes Analyzed	Alphanumeric		48018-6
617167	Disclaimer	Alphanumeric		62364-5
617168	Released By	Alphanumeric		18771-6

**LOINC®** and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

Supplemental

### **CPT Code Information:**

81410

### **Reference Values:**

An interpretive report will be provided.