
Reporting Title: Congenital Heart Disease Gene Panel

Performing Location: Rochester

Ordering Guidance:

Chromosomal microarray is often used as a first-tier test in the setting of congenital heart disease. If chromosomal microarray testing is desired, order either CMACB / Chromosomal Microarray, Congenital, Blood or CMAP / Chromosomal Microarray, Prenatal, Amniotic Fluid/Chorionic Villus Sampling.

Customization of this panel and single gene analysis for any gene present on this panel are available. For more information see CGPH / Custom Gene Panel, Hereditary, Next-Generation Sequencing, Varies.

Targeted testing for familial variants (also called site-specific or known mutations testing) is available for the genes on this panel. See FMTT / Familial Variant, Targeted Testing, Varies. To obtain more information about this testing option, call 800-533-1710.

Additional Testing Requirements:

All prenatal specimens must be accompanied by a maternal blood specimen; order MATCC / Maternal Cell Contamination, Molecular Analysis, Varies on the maternal specimen as this must be a different order number than the prenatal specimen.

Shipping Instructions:

Specimen preferred to arrive within 96 hours of collection.

Necessary Information:

[Prior Authorization](#) is available, **but not required**, for this test. If proceeding with the prior authorization process, submit the required form with the specimen.

Specimen Requirements:

Patient Preparation: A previous bone marrow transplant from an allogenic donor will interfere with testing. For instructions for testing patients who have received a bone marrow transplant, call 800-533-1710.

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Any anticoagulant

Specimen Volume: 3 mL

Collection Instructions:

1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. **Do not aliquot.**

Specimen Stability Information: Ambient (preferred)/Refrigerated

Due to its complexity, consultation with the laboratory is required for all prenatal testing; call 800-533-1710 to speak to a genetic counselor.

Specimen Type: Amniotic fluid
Container/Tube: Amniotic fluid container
Specimen Volume: 20 mL
Specimen Stability Information: Refrigerated (preferred)/Ambient
Additional information:
1. If amniotic fluid or nonconfluent cultures are received, CULAF / Culture for Genetic Testing, Amniotic Fluid will be added at an additional charge.
2. **All prenatal specimens must be accompanied by a maternal blood specimen;** order MATCC / Maternal Cell Contamination, Molecular Analysis, Varies on the maternal specimen.

Specimen Type: Chorionic villi
Container/Tube: 15-mL tube containing 15 mL of transport media
Specimen Volume: 20 mg
Specimen Stability Information: Refrigerated
Additional Information:
1. If nonconfluent cultures are received, CULFB / Fibroblast Culture for Biochemical or Molecular Testing will be added at an additional charge.
2. **All prenatal specimens must be accompanied by a maternal blood specimen;** order MATCC / Maternal Cell Contamination, Molecular Analysis, Varies on the maternal specimen.

Acceptable:
Specimen Type: Confluent cultured cells
Container/Tube: T-25 flask
Specimen Volume: 2 Flasks
Collection Instructions: Submit confluent cultured cells from another laboratory.
Specimen Stability Information: Ambient (preferred)/Refrigerated
Additional Information: **All prenatal specimens must be accompanied by a maternal blood specimen;** order MATCC / Maternal Cell Contamination, Molecular Analysis, Varies on the maternal specimen.

Forms:
1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available:
-[Informed Consent for Genetic Testing](#) (T576)
-[Informed Consent for Genetic Testing \(Spanish\)](#) (T826)
2. [Congenital Heart Disease Genetic Testing Patient Information](#)
3. [Congenital Heart Disease Gene Panel \(CHDGG\) Prior Authorization Ordering Instructions](#)
4. If not ordering electronically, complete, print, and send a [Cardiovascular Test Request](#) (T724) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
-----------	----------------	------	------	--------

617198	Test Description	Alphanumeric		62364-5
617199	Specimen	Alphanumeric		31208-2
617200	Source	Alphanumeric		31208-2
617201	Result Summary	Alphanumeric		50397-9
617202	Result	Alphanumeric		82939-0
617203	Interpretation	Alphanumeric		69047-9
617204	Additional Results	Alphanumeric		82939-0
617205	Resources	Alphanumeric		99622-3
617206	Additional Information	Alphanumeric		48767-8
617207	Method	Alphanumeric		85069-3
617208	Genes Analyzed	Alphanumeric		48018-6
617209	Disclaimer	Alphanumeric		62364-5
617210	Released By	Alphanumeric		18771-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

Supplemental

CPT Code Information:

- 81404
- 81405 x 3
- 81406 x 6
- 81407 x 3
- 81408
- 81479
- 81479 (if appropriate for government payers)
- 81265-Maternal cell contamination (if appropriate)
- 88233-Tissue culture, skin, solid tissue biopsy (if appropriate)
- 88235-Amniotic Fluid culture (if appropriate)
- 88240-Cryopreservation (if appropriate)

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
CULFB	Fibroblast Culture for Genetic Test	1	88233	No	Yes
CULAF	Amniotic Fluid Culture/Genetic Test	1	88235	No	Yes
MATCC	Maternal Cell Contamination, B	1	81265	No	Yes
_STR1	Comp Analysis using STR (Bill only)	1	81265	No	No, (Bill only)
_STR2	Add'l comp analysis w/STR (Bill Only)	1	81266	No	No, (Bill only)

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
CULAF	52304	Result Summary	Alphanumeric		50397-9
CULAF	52306	Interpretation	Alphanumeric		69965-2
CULAF	52305	Result	Alphanumeric		82939-0
CULAF	CG767	Reason for Referral	Alphanumeric		42349-1
CULAF	52307	Specimen	Alphanumeric		31208-2
CULAF	52308	Source	Alphanumeric		31208-2
CULAF	52309	Method	Alphanumeric		85069-3
CULAF	54641	Additional Information	Alphanumeric		48767-8
CULAF	52310	Released By	Alphanumeric		18771-6
CULFB	52327	Result Summary	Alphanumeric		50397-9
CULFB	52329	Interpretation	Alphanumeric		69965-2
CULFB	52328	Result	Alphanumeric		82939-0
CULFB	CG770	Reason for Referral	Alphanumeric		42349-1
CULFB	CG899	Specimen	Alphanumeric		31208-2
CULFB	52331	Source	Alphanumeric		31208-2
CULFB	52332	Method	Alphanumeric		85069-3
CULFB	54625	Additional Information	Alphanumeric		48767-8
CULFB	52333	Released By	Alphanumeric		18771-6
MATCC	53285	Result Summary	Alphanumeric		50397-9
MATCC	53286	Result	Alphanumeric		40704-9
MATCC	53287	Interpretation	Alphanumeric		69047-9
MATCC	53288	Reason for referral	Alphanumeric		42349-1
MATCC	53289	Specimen	Alphanumeric		31208-2
MATCC	53290	Source	Alphanumeric		31208-2
MATCC	53291	Released By	Alphanumeric		18771-6
MATCC	55150	Method	Alphanumeric		85069-3

Reference Values:

An interpretive report will be provided