

Reporting Title: Orthostatic Proteinuria, Random, U

Performing Location: Rochester

Specimen Requirements:

Patient Preparation: Specimens should be collected before fluorescein is given or not collected until at least 24 hour later.

Container/Tube: Plastic, 5-mL tube

Specimen Volume: 4 mL

Collection Instructions:

1. This is a nighttime (supine) collection.
2. Void immediately prior to lying down to sleep.
3. Collect an 8-hour (nighttime) urine specimen.
4. Upon awaking in the morning, the patient should immediately void into the specimen cup.
4. No preservative.
5. Invert well before taking 4-mL (maximum) aliquot.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	30 days	
	Ambient	24 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
CRETR	Creatinine, Random, U	Numeric	mg/dL	2161-8
PRCON	Protein, Total, Random, U	Numeric	mg/dL	2888-6
RATO2	Protein/Creatinine Ratio	Numeric	mg/mg	2890-2

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
PRCON	Protein, Total, Random, U	1	84156	Yes	No
RATO2	Protein/Creatinine Ratio			Yes	No
CRETR	Creatinine, Random, U	1	82570	Yes	No

CPT Code Information:

CRETR- 82570

PRCON- 84156

Reference Values:

PROTEIN/CREATININE RATIO:

> or =18 years: <0.18 mg/mg creatinine

CREATININE:> or =18 years = 16-326 mg/dL

Reference values have not been established for patients younger than 18 years of age.