

---

**Reporting Title:** Multiple sulfatase deficiency, BS**Performing Location:** Rochester**Necessary Information:**

1. Patient's age is required.
2. Reason for testing is required.

**Specimen Requirements:****Submit only 1 of the following specimen types:****Preferred:****Specimen Type:** Blood spot**Supplies:** Card-Blood Spot Collection (Filter Paper) (T493)**Container/Tube:****Preferred:** Blood Spot Collection Card**Acceptable:** Whatman Protein Saver 903 Paper, PerkinElmer 226 filter paper, Munktell filter paper, or blood collected in tubes containing ACD or EDTA and dried on filter paper.**Specimen Volume:** 2 Blood spots**Collection Instructions:**

1. An alternative blood collection option for a patient 1 year of age or older is a fingerstick. For detailed instructions, see [How to Collect Dried Blood Spot Samples](#).
2. At least 2 spots should be complete (ie, unpunched).
3. Let blood dry on filter paper at room temperature in a horizontal position for a minimum of 3 hours.
4. Do not expose specimen to heat or direct sunlight.
5. Do not stack wet specimens.
6. Keep specimen dry.

**Specimen Stability Information:** Refrigerated (preferred) 60 days/Ambient 7 days/Frozen 60 days**Additional Information:**

1. For collection instructions, see [Blood Spot Collection Instructions](#)
2. For collection instructions in Spanish, see [Blood Spot Collection Card-Spanish Instructions](#) (T777)
3. For collection instructions in Chinese, see [Blood Spot Collection Card-Chinese Instructions](#) (T800)

**Acceptable:****Specimen Type:** Whole Blood**Container/Tube:****Preferred:** Lavender top (EDTA)**Acceptable:** Yellow top (ACD)**Specimen Volume:** 2 mL**Collection Instructions:** Send whole blood specimen in original tube. **Do not aliquot.****Specimen Stability Information:** Refrigerate (preferred) 7 days/Ambient 48 hours**Forms:**

1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available:

[-Informed Consent for Genetic Testing](#) (T576)[-Informed Consent for Genetic Testing-Spanish](#) (T826)

2. [Biochemical Genetics Patient Information](#) (T602)
3. If not ordering electronically, complete, print, and send a [Biochemical Genetics Test Request](#) (T798) with the specimen.

| Specimen Type | Temperature | Time | Special Container |
|---------------|-------------|------|-------------------|
| Whole blood   | Varies      |      |                   |

Ask at Order Entry (AOE) Questions:

| Test ID | Question ID | Description                                                                                                                                                                         | Type        | Reportable |
|---------|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|------------|
| MSDBS   | BG755       | Reason for Referral: <ul style="list-style-type: none"><li>• Rule out MSD</li><li>• Follow up of known MSD</li><li>• Multiple Sulfatase deficiency</li><li>• Not Provided</li></ul> | Answer List | Yes        |

Result Codes:

| Result ID | Reporting Name                    | Type         | Unit      | LOINC®   |
|-----------|-----------------------------------|--------------|-----------|----------|
| BG755     | Reason for Referral               | Alphanumeric |           | 42349-1  |
| 618430    | Iduronate-2-sulfatase             | Numeric      | nmol/mL/h | 79462-8  |
| 618431    | Heparan-N-sulfatase               | Numeric      | nmol/mL/h | 104113-6 |
| 618432    | N-acetylgalactosamine-6-sulfatase | Numeric      | nmol/mL/h | 88019-5  |
| 618433    | Interpretation                    | Alphanumeric |           | 59462-2  |
| 618429    | Reviewed By                       | Alphanumeric |           | 18771-6  |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82657  
83864 (if appropriate)

Reflex Tests:

| Test Id | Reporting Name            | CPT Units | CPT Code | Always Performed | Available Separately |
|---------|---------------------------|-----------|----------|------------------|----------------------|
| MPSBS   | Mucopolysaccharidosis, BS | 1         | 83864    | No               | Yes                  |

Result Codes for Reflex Tests:

| Test ID | Result ID | Reporting Name | Type | Unit | LOINC® |
|---------|-----------|----------------|------|------|--------|
|---------|-----------|----------------|------|------|--------|

|       |        |                        |              |        |         |
|-------|--------|------------------------|--------------|--------|---------|
| MPSBS | 43693  | Dermatan Sulfate       | Numeric      | nmol/L | 90233-8 |
| MPSBS | 43694  | Heparan Sulfate        | Numeric      | nmol/L | 90235-3 |
| MPSBS | 43695  | Interpretation (MPSBS) | Alphanumeric |        | 59462-2 |
| MPSBS | 43696  | Reviewed By            | Alphanumeric |        | 18771-6 |
| MPSBS | BA2869 | Total Keratan Sulfate  | Alphanumeric | nmol/L | 90236-1 |

Reference Values:

Iduronate-2-sulfatase: >4.30 nmol/mL/hour  
Heparan-N-sulfatase: >0.06 nmol/mL/hour  
N-acetylgalactosamine-6-sulfatase: >0.70 nmol/mL/hour

An interpretive report will be provided.