

# **Test Definition: BPRPV**

Bordetella pertussis and Bordetella parapertussis, Molecular Detection, PCR, Varies

Reporting Title: Bordetella, PCR, Varies

Performing Location: Rochester

**Necessary Information:** Specimen source is required.

## **Specimen Requirements:**

The high sensitivity of amplification by polymerase chain reaction (PCR) requires the specimen to be processed in an environment in which contamination of the specimen by *Bordetella pertussis* or *Bordetella parapertussis* DNA is unlikely.

#### Submit only 1 of the following specimens:

**Preferred:** 

Specimen Type: Nasopharyngeal swab

**Supplies:** 

Culture Swab - Liquid Stuarts/Single Swab (NP Swab) (T515)

**Container/Tube:** Rayon swab with an aluminum shaft placed in transport medium such as a green-top nasopharyngeal swab (rayon mini-tip) with Stuart's media (no charcoal), or Stuart's media with charcoal, or Amies media with or without charcoal (Transwab Nasopharyngeal with Charcoal System).

#### **Additional Information:**

- 1. Swab transport containers without charcoal must contain a pledget saturated with either Stuart's or Amies liquid media. Clear semi-solid/solid media is gel and will be rejected.
- 2. Other swab or media types may be inhibitory to PCR testing and will be rejected.

#### Acceptable:

**Specimen Type:** Nasopharyngeal (not throat) aspirate/wash or nasal aspirate/wash **Container/Tube:** Sterile container with a screw top cap (no transport media)

Specimen Volume: Entire collection

#### Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T244) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Ambient	7 days	
	Frozen	7 days	

#### Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
BPRPV	BPRS	Specimen source:	Answer List	Yes
		<ul> <li>Nasopharyngeal swab</li> </ul>		



# **Test Definition: BPRPV**

Bordetella pertussis and Bordetella parapertussis, Molecular Detection, PCR, Varies

ı		_	
		<ul> <li>Nasal aspirate/wash</li> </ul>	
		<ul> <li>Nasopharyngeal aspirate/wash</li> </ul>	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
BPRS	Specimen source	Alphanumeric		31208-2
618312	Bordetella pertussis PCR	Alphanumeric		43913-3
618313	Bordetella parapertussis PCR	Alphanumeric		42588-4

LOINC® and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

#### **CPT Code Information:**

87798 x 2

87999 (if appropriate for government payers)

#### **Reference Values:**

Not applicable