

---

**Reporting Title:** PROS1 Gene, Full Gene NGS**Performing Location:** Rochester**Ordering Guidance:**

This test should only be considered if clinical and family history, initial coagulation screens, and/or initial antigen and activity tests suggest a diagnosis of protein S deficiency (see Testing Algorithm).

This test does not measure protein S activity or antigen levels.

-For assessment of free protein S activity, order S\_FX / Protein S Activity, Plasma.

-For assessment of plasma free protein S antigen, order PSTF / Protein S Antigen, Plasma.

For patients in whom hereditary protein S deficiency is strongly suspected and the plasma free protein S antigen level is normal, consider testing free protein S activity for detecting type II protein S deficiency, which is very rare. Order S\_FX / Protein S Activity, Plasma.

If genetic testing for hereditary blood clotting disorders using a larger panel is desired, a 16-gene comprehensive thrombosis panel is available; order GNTHR / Thrombosis Disorders, Comprehensive Gene Panel, Next-Generation Sequencing, Varies.

Testing for the *PROS1* gene as part of a customized panel is available. For more information see CGPH / Custom Gene Panel, Hereditary, Next-Generation Sequencing, Varies.

Targeted testing for familial variants (also called site-specific or known variants testing) is available for the *PROS1* gene. See FMTT / Familial Variant, Targeted Testing, Varies. To obtain more information about this testing option, call 800-533-1710.

**Shipping Instructions:**

Specimen preferred to arrive within 96 hours of collection.

**Necessary Information:**

[Rare Coagulation Disorder Patient Information](#) is required. Testing may proceed without the patient information, however, the information aids in providing a more thorough interpretation. Ordering providers are strongly encouraged to fill out the form and send with the specimen.

**Specimen Requirements:**

**Specimen Type:** Whole blood

**Patient Preparation:** A previous bone marrow transplant from an allogenic donor will interfere with testing. Call 800-533-1710 for instructions for testing patients who have received a bone marrow transplant.

**Container/Tube:**

**Preferred:** Lavender top (EDTA)

**Acceptable:** Yellow top (ACD)

**Specimen Volume:** 3 mL

**Collection Instructions:**

1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. **Do not aliquot.**

# Test Definition: GNPRS

Protein S Deficiency, PROS1 Gene,  
Next-Generation Sequencing, Varies

**Specimen Stability Information:** Ambient (preferred) 4 days/Refrigerated

**Forms:**

- 1. [Rare Coagulation Disorder Patient Information \(T824\)](#) is **required**.
- 2. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available:
  - [Informed Consent for Genetic Testing \(T576\)](#)
  - [Informed Consent for Genetic Testing \(Spanish\) \(T826\)](#)
- 3. [If not ordering electronically, complete, print, and send an Coagulation Test Request \(T753\)](#) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
619188	Test Description	Alphanumeric		62364-5
619189	Specimen	Alphanumeric		31208-2
619190	Source	Alphanumeric		31208-2
619191	Result Summary	Alphanumeric		50397-9
619192	Result	Alphanumeric		82939-0
619193	Interpretation	Alphanumeric		59465-5
619194	Additional Results	Alphanumeric		82939-0
619195	Resources	Alphanumeric		99622-3
619196	Additional Information	Alphanumeric		48767-8
619197	Method	Alphanumeric		85069-3
619198	Genes Analyzed	Alphanumeric		82939-0
619199	Disclaimer	Alphanumeric		62364-5
619200	Released By	Alphanumeric		18771-6

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

Supplemental

**CPT Code Information:**

81479

**Reference Values:**

An interpretive report will be provided.