

# **Test Definition: HYDCU**

Hydrocodone with Metabolite Confirmation, Random, Urine

Reporting Title: Hydrocodone w/metabolite Conf, U

Performing Location: Rochester

#### **Ordering Guidance:**

For situations where chain of custody is required, a Chain of Custody Kit (T282) is available. For chain-of-custody testing, order OPATX / Opiates Confirmation, Chain of Custody, Random, Urine.

Additional drug panels and specific requests are available; call 800-533-1710 or 507-266-5700.

#### **Additional Testing Requirements:**

**If urine creatinine is required** or adulteration of the sample is suspected, order ADULT / Adulterants Survey, Random, Urine.

#### **Specimen Requirements:**

**Supplies:** Sarstedt 5 mL Aliquot Tube (T914) **Collection Container/Tube:** Plastic urine container **Submission Container/Tube:** Plastic, 5-mL tube

**Specimen Volume:** 3 mL Collection Instructions:

- 1. Collect a random urine specimen.
- 2. No preservative.

# **Additional Information:**

- 1. No specimen substitutions.
- 2. STATS are **not** accepted for this test.
- 3. Submitting <1 mL will compromise our ability to perform all necessary testing.

#### Forms:

If not ordering electronically, complete, print, and send a <u>Therapeutics Test Request</u> (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container	
Urine	Refrigerated (preferred)	28 days		
	Frozen	28 days		
	Ambient	14 days		

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
62614	Hydrocodone-by LC-MS/MS	Alphanumeric	ng/mL	16252-9
35966	Hydromorphone-by LC-MS/MS	Alphanumeric	ng/mL	16998-7
36026	Hydrocodone Interpretation	Alphanumeric		69050-3
41999	Norhydrocodone-by LC-MS/MS	Alphanumeric	ng/mL	61422-2

LOINC® and CPT codes are provided by the performing laboratory.



# **Test Definition: HYDCU**

Hydrocodone with Metabolite Confirmation, Random, Urine

# **Supplemental Report:**

No

#### **CPT Code Information:**

80361

G0480 (if appropriate)

# **Reference Values:**

Negative

Cutoff concentrations-liquid chromatography-tandem mass spectrometry:

Hydrocodone: 25 ng/mL Norhydrocodone: 25 ng/mL Hydromorphone: 25 ng/mL