

## **Test Definition: MPNR**

Myeloproliferative Neoplasm, JAK2 V617F with Reflex to CALR and MPL, Varies

Reporting Title: MPN (JAK2 V617F, CALR, MPL) Reflex

Performing Location: Rochester

#### **Shipping Instructions:**

Specimen must arrive within 7 days of collection.

#### **Necessary Information:**

## The following information is required:

- 1. Pertinent clinical history
- 2. Clinical or morphologic suspicion
- 3. Date of collection
- 4. Specimen source

## **Specimen Requirements:**

Submit only 1 of the following specimens:

Specimen Type: Blood

Container/Tube: Lavender top (EDTA) or yellow top (ACD solution B)

**Specimen Volume:** 3 mL **Collection Instructions:** 

1. Invert several times to mix blood.

2. Send specimen in original tube. **Do not** aliquot.

3. Label specimen as blood.

Specimen Stability Information: Ambient (preferred)/Refrigerate

**Specimen Type:** Bone marrow

Container/Tube: Lavender top (EDTA) or yellow top (ACD solution B)

**Specimen Volume:** 2 mL **Collection Instructions:** 

1. Invert several times to mix bone marrow.

2. Send specimen in original tube. **Do not** aliquot.

3. Label specimen as bone marrow.

Specimen Stability Information: Ambient (preferred)/Refrigerate

Specimen Type: Extracted DNA from blood or bone marrow

Container/Tube: 1.5- to 2-mL tube
Specimen Volume: Entire specimen

Collection Instructions: Label specimen as extracted DNA from blood or bone marrow and provide indication of volume

and concentration of the DNA.

Specimen Stability Information: Frozen (preferred)/Refrigerate/Ambient

## Forms:

If not ordering electronically, complete, print, and send a <u>Hematopathology/Cytogenetics Test Request</u> (T726) with the specimen.



# **Test Definition: MPNR**

Myeloproliferative Neoplasm, JAK2 V617F with Reflex to CALR and MPL, Varies

Specimen Type	Temperature	Time	Special Container
Varies	Varies	7 days	

## Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
MPNR	MP023	Specimen Type	Plain Text	Yes

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
36687	Final Diagnosis	Alphanumeric		22637-3
39725	MPNR Result	Alphanumeric		No LOINC Needed

LOINC® and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

## **CPT Code Information:**

81270-JAK2 (Janus kinase 2) (eg, myeloproliferative disorder) gene analysis, p.Val617Phe (V617F) variant 81219-CALR (calreticulin) (eg, myeloproliferative disorders), gene analysis, common variants in exon 9 (if appropriate) 81339-MPL (MPL proto-oncogene, thrombopoietin receptor) (eg, myeloproliferative disorder) gene analysis; sequence analysis, exon 10 (if appropriate)

## **Reflex Tests:**

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
CALX	CALR, Gene Mutation, Exon 9, Reflex	1	81219	No	No, (bill only)
MPLR	MPL Exon 10 Mutation Detection, R	1	81339	No	No, (bill only)

#### **Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
MPLR	36672	Final Diagnosis	Alphanumeric		22637-3
CALX	36998	Final Diagnosis	Alphanumeric		22637-3

#### **Reference Values:**

An interpretive report will be provided.