

Reporting Title: Fentanyl Screen w/Reflex, U

Performing Location: Rochester

Ordering Guidance:

For situations where chain of custody is required, a Chain of Custody Kit (T282) is available. For chain-of-custody testing, order FENTX / Fentanyl with Metabolite Confirmation, Chain of Custody, Urine.

For monitoring therapeutic drug levels, order FENTS / Fentanyl, Serum.

Additional Testing Requirements:

If urine creatinine is required or adulteration of the sample is suspected, the ADULT / Adulterants Survey, Random, Urine test should also be ordered.

Specimen Requirements:

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Clean, plastic urine collection container

Submission Container/Tube: Plastic, 5-mL tube

Specimen Volume: 5 mL

Collection Instructions:

1. Collect a random urine specimen.
2. No preservative.

Forms:

If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
63060	Fentanyl Screen, U	Alphanumeric	ng/mL	59673-4

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
FENS	Fentanyl Screen, U	1	80307	Yes	Yes

CPT Code Information:

80307

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
FENTU	Fentanyl w/metabolite Conf, U	1	80354	No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
FENTU	45531	Norfentanyl by LC-MS/MS	Alphanumeric	ng/mL	58383-1
FENTU	45530	Fentanyl by LC-MS/MS	Alphanumeric	ng/mL	58381-5
FENTU	45532	Fentanyl Interpretation	Alphanumeric		69050-3

Reference Values:

Negative

Screening cutoff concentration: 2 ng/mL