

Reporting Title: ETV6 (12p13.2), FISH, Ts
Performing Location: Rochester

Ordering Guidance:
This test does not include a pathology consultation. If a pathology consultation is requested, order PATHC / Pathology Consultation, and appropriate testing will be added at the discretion of the pathologist and performed at an additional charge.

Multiple oncology (cancer) gene panels are also available. For more information see [Hematology, Oncology, and Hereditary Test Selection Guide](#)

Additional Testing Requirements:
Confirmation testing by next-generation sequencing to resolve atypical or unbalanced fluorescence in situ hybridization results of this gene region is available, order SARCP / Sarcoma Targeted Gene Fusion/Rearrangement Panel, Next-Generation Sequencing, Tumor.

Shipping Instructions:
Advise Express Mail or equivalent if not on courier service.

Necessary Information:

1. **A pathology report is required for testing to be performed.** If not provided, appropriate testing and/or interpretation may be compromised or delayed. Acceptable pathology reports include working drafts, preliminary pathology, or surgical pathology reports.

2. **The following information must be included in the report provided.**

- Patient name
- Block number - **must be on all blocks, slides, and paperwork**
- Date of collection
- Tissue source

3. **A reason for testing must be provided.** If this information is not provided, an appropriate indication for testing may be entered by Mayo Clinic Laboratories.

Specimen Requirements:
Submit only 1 of the following specimens:

Preferred
Specimen Type: Tissue block
Collection Instructions: Submit a formalin-fixed, paraffin-embedded tumor tissue block. Blocks prepared with alternative fixation methods may be acceptable; provide fixation method used.

Acceptable
Specimen Type: Tissue slides
Slides: 1 Hematoxylin and eosin stained and 4 unstained
Collection Instructions: Submit 1 slide stained with hematoxylin and eosin and 4 consecutive unstained, positively-charged, unbaked slides with 5 micron-thick sections of the tumor tissue.

Forms:
If not ordering electronically, complete, print, and send an [Oncology Test Request](#) (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container
---------------	-------------	------	-------------------

Tissue	Ambient (preferred)		
	Refrigerated		

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
ETV6F	CG977	Reason for Referral	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
46837	Result Summary	Alphanumeric		50397-9
46838	Interpretation	Alphanumeric		69047-9
46839	Result	Alphanumeric		62356-1
46840	Reason For Referral	Alphanumeric		42349-1
46841	Specimen	Alphanumeric		31208-2
46842	Source	Alphanumeric		39111-0
46843	Tissue ID	Alphanumeric		80398-1
46844	Method	Alphanumeric		85069-3
46845	Additional Information	Alphanumeric		48767-8
46846	Disclaimer	Alphanumeric		62364-5
46847	Released By	Alphanumeric		18771-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

- 88271x2
- 88291 -DNA probe, each (first probe set), Interpretation and report
- 88271x2 -DNA probe, each; each additional probe set (if appropriate)
- 88271x1 -DNA probe, each; coverage for sets containing 3 probes (if appropriate)
- 88271x2 -DNA probe, each; coverage for sets containing 4 probes (if appropriate)
- 88271x3 -DNA probe, each; coverage for sets containing 5 probes (if appropriate)
- 88274 w/modifier 52 -Interphase in situ hybridization, <25 cells, each probe set (if appropriate)
- 88274 -Interphase in situ hybridization, 25 to 99 cells, each probe set (if appropriate)
- 88275 -Interphase in situ hybridization, 100 to 300 cells, each probe set (if appropriate)

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
_PBCT	Probe, +2	1	88271	No	No

_PADD	Probe, +1	1	88271	No	No
_PB02	Probe, +2	1	88271	No	No
_PB03	Probe, +3	1	88271	No	No
_IL25	Interphases, <25	1	88274	No	No
_I099	Interphases, 25-99	1	88274	No	No
_I300	Interphases, >=100	1	88275	No	No

Reference Values:

An interpretive report will be provided.