

Reporting Title: IgA Heavy Light Chains (HLC), S
Performing Location: LabCorp Burlington

Specimen Requirements:

Specimen Type: Serum
Container/Tube: Red top or SST
Specimen Volume: 0.75 mL

Collection Instructions: Draw blood in a plain red-top tube, serum gel tube(s) is also acceptable. Separate serum immediately after coagulation (30 minutes) to prevent hemolysis. Send 0.75 mL of serum refrigerated in a plastic vial.

NOTE: Patient should be fasting for eight hours to avoid lipemic sample interference.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Ambient	14 days	
	Frozen	14 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z5616	IgA Kappa, S	Alphanumeric		74864-0
Z5617	IgA Lambda, S	Alphanumeric		74865-7
Z5618	IgA K/L HLC Ratio	Alphanumeric		74869-9

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83883 x 2

Reference Values:

IgA Kappa (g/L): 0.48-2.82
IgA Lambda (g/L): 0.36-1.98
IgA Kappa:IgA Lambda ratio: 0.80-2.04