

Test Definition: FNAD

Nadolol, Serum/Plasma

Reporting Title: Nadolol
Performing Location: NMS Labs

Specimen Requirements:

Must submit one specimen per order. Specimens cannot be shared between multiple orders.
Submit only 1 of the following specimens.

Serum

Specimen Type: Serum

Collection Container/Tube: Red top

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain, red-top tube(s). (Serum gel tube is not acceptable.) Spin down and send 1

mL of serum refrigerated in a plastic preservative-free vial.

Plasma

Specimen Type: Plasma

Collection Container/Tube: Lavender top or pink top (EDTA)

Specimen volume: 1 mL

Collection Instructions: Draw blood in a lavender-top (EDTA) or Pink top tube(s). (Plasma gel tube is not acceptable.)

Spin down and send 1 mL of EDTA plasma refrigerated in a plastic preservative-free vial.

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	30 days	
	Ambient	30 days	
	Frozen	30 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z1430	Nadolol	Alphanumeric	ng/mL	12406-5
Z1867	Reporting Limit	Alphanumeric	ng/mL	19147-8

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80299

Reference Values:

Reporting limit determined each analysis



Test Definition: FNAD

Nadolol, Serum/Plasma

None Detected ng/mL

Mean steady-state plasma levels following a daily regimen:

80 mg: 26 - 36 ng/mL 160 mg: 52 - 74 ng/mL 320 mg: 154 - 191 ng mL