

Reporting Title: Inhibin A, Tumor Marker, S

Performing Location: Rochester

**Ordering Guidance:**

For the initial evaluation of patients suspected of having a granulosa cell tumor of the ovary, order INHAB / Inhibin A and B, Tumor Marker, Serum. If the results of the profile show that either inhibin A or B are elevated, consider monitoring the patient with the individual tests, INHA / Inhibin A, Tumor Marker, Serum or INHB / Inhibin B, Serum.

**Specimen Requirements:**

**Supplies:** Sarstedt Aliquot Tube 5 mL (T914)

**Collection Container/Tube:**

**Preferred:** Serum gel

**Acceptable:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 0.6 mL

**Collection Instructions:** Centrifuge and aliquot serum into a plastic vial.

**Forms:**

[If not ordering electronically, complete, print, and send an Oncology Test Request](#) (T729) with the specimen.

| Specimen Type | Temperature              | Time    | Special Container |
|---------------|--------------------------|---------|-------------------|
| Serum         | Refrigerated (preferred) | 14 days |                   |
|               | Frozen                   | 90 days |                   |
|               | Ambient                  | 7 days  |                   |

**Result Codes:**

| Result ID | Reporting Name             | Type    | Unit  | LOINC®  |
|-----------|----------------------------|---------|-------|---------|
| INHA      | Inhibin A, Tumor Marker, S | Numeric | pg/mL | 23883-2 |

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86336

**Reference Values:**

Males: <5.0 pg/mL

Females

<11 years: <5.0 pg/mL

11-17 years: <98 pg/mL

---

Premenopausal: <98 pg/mL  
Postmenopausal: <5.0 pg/mL