

Test Definition: PIPA

Pipecolic Acid, Serum

Reporting Title: Pipecolic Acid, S **Performing Location:** Rochester

Necessary Information:

Patient's age is required.

Specimen Requirements:

Patient Preparation: Fasting 12 hours or more. (Collect specimens from infants and small children just before next feeding)
Supplies: Sarstedt Aliquot Tube, 5 mL (T914)
Collection Container/Tube:
Preferred: Serum gel
Acceptable: Red top
Submission Container/Tube: Plastic vial
Specimen Volume: 1.5 mL
Collection Instructions: Centrifuge and aliquot serum into plastic vial.

Forms:

If not ordering electronically, complete, print, and send a Biochemical Genetics Test Request (T798) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	94 days	
	Refrigerated	14 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
81326	Pipecolic Acid, S	Numeric	nmol/mL	32334-5
29964	Reviewed By	Alphanumeric		18771-6
29962	Interpretation	Alphanumeric		59462-2

LOINC[®] and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82542

Reference Values:

<6 months: < or =6.0 nmol/mL 6 months-<1 year: < or =5.9 nmol/mL 1-17 years: < or =4.3 nmol/mL > or =18 years: < or =7.4 nmol/mL



Test Definition: PIPA

Pipecolic Acid, Serum