

Test Definition: FS

Fungal Smear, Varies

Reporting Title: Fungal Smear **Performing Location:** Rochester

Shipping Instructions:

Specimen should arrive within 24 hours of collection.

Necessary Information:

Specimen source is required.

Specimen Requirements:

Preferred Specimen Type: Body fluid Container/Tube: Sterile container Specimen Volume: Entire collection

Preferred Specimen Type: Fresh tissue Container/Tube: Sterile container Specimen Volume: Pea size

Collection Instructions: Tissue should be placed in small amount of sterile saline or sterile water.

Specimen Type: Bone marrow

Container/Tube: Sterile container, SPS/Isolator system, or green top (lithium or sodium heparin)

Specimen Volume: Entire collection

Specimen Type: Respiratory specimen Container/Tube: Sterile container Specimen Volume: Entire collection

Specimen Type: Urine

Container/Tube: Sterile container

Specimen Volume: 2 mL

Collection Instructions: Collect a random urine specimen.

Acceptable Specimen Type: Swab

Sources: Dermal, ear, mouth, ocular, throat, or wound

Container/Tube: Culture transport swab (non-charcoal) Culturette or Eswab

Specimen Volume: Swab Collection Instructions:

- 1. Before collecting specimen, wipe away any excessive amount of secretion and discharge, if appropriate.
- 2. Obtain secretions or fluid from source with sterile swab.
- 3. If smear and culture are requested or both a bacterial culture and fungal culture are requested, collect a second swab to maximize test sensitivity.

Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T244) with the specimen.



Test Definition: FS

Fungal Smear, Varies

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Ambient	7 days	

Ask at Order Entry (AOE) Questions:

Test ID Question ID		Description	Туре	Reportable
FS	Q00M0023	Specimen Source	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
FS	Fungal Smear	Alphanumeric		658-5

 $\ensuremath{\mathsf{LOINC}} \ensuremath{^{\circledR}}$ and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

872N6

87176-Tissue processing (if appropriate)

Reflex Tests:

	Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
ĺ	TISSR	Tissue Processing	1	87176	No	No, (Bill Only)

Reference Values:

Negative