

Reporting Title: Sertraline (Zoloft)
Performing Location: Medtox Laboratories, Inc.

Specimen Requirements:
Submit only 1 of the following specimens:

Plasma
Specimen Type: Plasma
Container/Tube: Green Top
Specimen Volume: 2 mL
Collection Instructions: Draw blood in a green-top (sodium heparin) tube, **plasma gel tube is not acceptable**. Spin down and send 2 mL of sodium heparin plasma refrigerated in a plastic vial.

Serum
Specimen Type: Serum
Container/Tube: Red
Specimen Volume: 2 mL
Collection Instructions: Draw blood in a plain, red-top tube, **serum gel tube is not acceptable**. Spin down and send 2 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Frozen	180 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z2298	Sertraline	Alphanumeric		6906-2
Z2299	Desmethylsertraline	Alphanumeric		6897-3

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
80332

Reference Values:
Sertraline:

Reference Range: 30 - 200 ng/mL
Report Limit 10 ng/mL

Desmethylsertraline: ng/mL
No reference range provided

The stated reference range is the range of observed steady-state concentrations in individuals receiving therapeutic dosage regimens of sertraline. This is not a defined therapeutic range.
Report Limit 10 ng/mL