

Reporting Title: Phenosense Entry Fuzeon
Performing Location: Monogram Biosciences, Inc

Specimen Requirements:

Specimen Type: Plasma
Container/Tube: Lavender-top (EDTA) or plasma preparation tube (PPT)
Specimen Volume: 3 mL
Collection Instructions: Draw blood into **two** 5-mL PPT (pearl top) or EDTA (lavender top) tube. Centrifuge specimen within six hours of collection. Transfer plasma to one or more polypropylene screw-capped tube(s) and freeze. Freeze immediately. Send 3 mL plasma in a screw-cap vial frozen.

To avoid delays in turnaround time when requesting multiple tests, please submit separate frozen specimens for each test requested.

RECOMMENDED:

- 1. Patient's most recent viral load
- 2. Viral load collection date

NOTE:

- 1. Intended for use only for patients with viral loads greater than or equal to 500 copies/mL. For best results, viral loads should be confirmed within two weeks prior to submission for testing at Monogram.
- 2. Patient samples submitted <30 days apart are considered duplicate and will be canceled.

Specimen Type	Temperature	Time	Special Container
Plasma EDTA	Frozen	14 days	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
FPFUZ	Z2098	Most Recent Viral Load	Plain Text	Yes
FPFUZ	Z2099	Viral Load Collected Date	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z1041	Phenosense Entry Fuzeon	Alphanumeric		45182-3

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

Referral

CPT Code Information:

87903

Reference Values:

A final report will be attached in MayoAccess.