

Barley IgG

Overview

#### Method Name

Enzyme Immunoassay (FEIA)

## NY State Available

No

## Specimen

Specimen Type Serum

### **Specimen Required**

Specimen Type: Serum

Container/Tube: Red Top

Specimen Volume: 0.5 mL

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

## Specimen Minimum Volume

0.5 mL

**Reject Due To** 

### **Specimen Stability Information**

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Refrigerated (preferred) | 28 days  |                   |
|               | Frozen                   | 365 days |                   |
|               | Ambient                  | 7 days   |                   |

# Clinical & Interpretive

### **Reference Values**

Reference ranges have not been established for food-specific IgG tests.

#### Interpretation

The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations



to select foods for evaluation by diet eliminations and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of this test.

# Performance

PDF Report No

Day(s) Performed Monday through Friday

**Report Available** 5 to 7 days

**Performing Laboratory Location** Eurofins Viracor

Fees & Codes

#### Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

### **Test Classification**

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

### **CPT Code Information**

86001

FBARG

### LOINC<sup>®</sup> Information

| Test ID   | Test Order Name  | Order LOINC <sup>®</sup> Value  |
|-----------|------------------|---------------------------------|
| FBARG     | Barley IgG       | 51917-3                         |
|           |                  |                                 |
| Result ID | Test Result Name | Result LOINC <sup>®</sup> Value |

51917-3

Barley IgG