

**Overview****Useful For**

Aiding in the diagnosis of *Kingella kingae* infection using whole blood specimens

**Method Name**

Real-Time Polymerase Chain Reaction (PCR)

**NY State Available**

Yes

**Specimen****Specimen Type**

Whole Blood EDTA

**Specimen Required**

The high sensitivity of amplification by PCR requires the specimen to be processed in an environment in which contamination of the specimen by *Kingella kingae* DNA is unlikely.

**Container/Tube:**

**Preferred:** Lavender top (EDTA)

**Acceptable:** Royal blue top (EDTA), pink top (EDTA), or sterile vial containing EDTA-derived aliquot

**Specimen Volume:** 1 mL

**Collection Instructions:** Send specimen in original tube (preferred).

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Whole Blood EDTA	Refrigerated (preferred)	7 days	
	Frozen	7 days	

**Clinical and Interpretive****Clinical Information**

*Kingella kingae* is a fastidious short Gram-negative bacillus that may colonize the oropharynx of young children. Colonization may occasionally lead to invasive disease via hematogenous dissemination, primarily in children younger than 4 years of age. This most commonly results in bone and joint infection; *K kingae* is the most frequent cause of osteomyelitis and septic arthritis in children aged 6 to 36 months. *K kingae* may also cause endocarditis, involving both native and prosthetic valves, in patients of any age and is considered part of the HACEK (*Haemophilus* species, *Aggregatibacter* species, *Cardiobacterium hominis*, *Eikenella corrodens*, and *Kingella* species) group of organisms, known for causing culture-negative endocarditis. *K kingae* produces a repeat-in-toxin (RTX) toxin.

Diagnosis of *K kingae* infection may be challenging due to the fastidious nature of the organism in culture. Evaluation of blood by PCR is a useful tool for the diagnosis of some cases of *K kingae* infection.

### Reference Values

Not applicable

### Interpretation

A positive result indicates the presence of *Kingella kingae* DNA.

A negative result indicates the absence of detectable *K kingae* DNA, but does not negate the presence of the organism and may occur due to inhibition of PCR, sequence variability underlying primers or probes, or the presence of *K kingae* DNA in quantities less than the limit of detection of the assay.

### Cautions

Test results should be used as an aid in diagnosis. The single assay should not be used as the only criteria to form a clinical conclusion, but results should be correlated with patient symptoms and clinical presentation. A negative result does not negate the presence of the organism or active disease.

This assay does not detect species of *Kingella* other than *kingae* or *negevensis* (see Supportive Data).

This assay cross-reacts with *Kingella negevensis*.<sup>(1)</sup>

### Supportive Data

This assay was validated by testing 30-spiked positive EDTA whole blood samples and 10-negative samples. No PCR inhibitors were encountered. The assay was 100% sensitive and specific. The assay showed no cross-reactivity when tested with a panel of 67 bacterial isolates, including *Kingella* species other than *kingae*. The limit of detection (LoD) in EDTA-whole blood was 1.3 CFU/mL.

### Clinical Reference

1. El Houmami N, Bzdreng J, Durand GA, et al: Molecular tests that target the RTX locus do not distinguish between *Kingella kingae* and the recently described *Kingella negevensis* species. *J Clin Microbiol* 2017;55:3113-3122
2. Murphy TF: In Mandell, Douglas, and Bennett's Principles and Practice of Infectious Diseases. Edited by GL Mandell, JE Bennett, R Dolin. Seventh edition. Philadelphia, Churchill Livingstone/Elsevier, 2010, pp 2774-2776
3. Zbinden R: *Aggregatibacter*, *Capnocytophaga*, *Eikenella*, *Kingella*, *Pasteurella*, and Other Fastidious or Rarely Encountered Gram-Negative Rods. In Manual of Clinical Microbiology. Edited by JH Jorgensen, KC Carroll, G Funke, MA Pfaller. 11th edition. Washington DC. ASM Press 2015, pp 652-666
4. Yagupsky P: *Kingella kingae*: carriage, transmission, and disease. [Clin Microbiol Rev](#). 2015 Jan;28(1):54-79
5. Madigan T, Cunningham SA, Ramanan P, et al: Real-Time PCR Assay for Detection of *Kingella kingae* in Children. *J Pediatr Infect Dis* 2018;13:216-233

## Performance

### Method Description

Nucleic acid is extracted from the specimen using the automated MagNA Pure instrument. Target specific primers are used to amplify the *txtB* gene region of *Kingella kingae*; amplification is monitored by detecting fluorescence produced by target specific FRET hybridization probes. This real-time PCR reaction takes place on a LightCycler instrument. Detection of the *K kingae* target is performed through melting curve analysis using the LightCycler software. (Cockerill FR, Uhl JR: Applications and challenges of real-time PCR for the clinical microbiology laboratory. In Rapid Cycle Real-Time PCR Methods and Applications. Edited by U Reischl, C Wittwer, F Cockerill. Berlin, Germany, Springer-Verlag, 2002, pp 3-27; Zbinden R: *Aggregatibacter*, *Capnocytophaga*, *Eikenella*, *Kingella*, *Pasteurella*, and Other Fastidious or Rarely Encountered Gram-Negative Rods. In Manual of Clinical Microbiology. 12th edition. Edited by K Carroll, M Pfaller. Washington DC, ASM Press, 2019, pp 656-669)

### PDF Report

No

### Day(s) Performed

Monday, Wednesday, Friday

### Report Available

2 to 7 days

### Specimen Retention Time

1 week

### Performing Laboratory Location

Rochester

## Fees and Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed, and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food and Drug Administration.

### CPT Code Information

87798

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
KKBRP	Kingella kingae PCR, B	65809-6

---

Result ID	Test Result Name	Result LOINC Value
48451	Specimen Source	31208-2
48338	Kingella kingae PCR, B	65809-6