

**Overview****Method Name**

ImmunoCAP FEIA

**NY State Available**

Yes

**Specimen****Specimen Type**

Serum

**Specimen Required**

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Gross hemolysis:	Mild OK; Gross Reject
Thawing:	Warm OK; Cold OK
Gross lipemia:	Reject
Gross icterus:	NA
Other:	NA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

**Clinical and Interpretive****Reference Values**

&lt;0.35 kU/L

**Interpretation**

Class	IgE (kU/L)	Comment
0	<0.35	Below Detection
1	0.35 - 0.69	Low Positive
2	0.70 - 3.49	Moderate Positive
3	3.50 - 17.49	Positive
4	17.50 - 49.99	Strong Positive
5	50.00 - 99.99	Very Strong Positive
6	>99.99	Very Strong Positive

**Performance****Method Description**

The ImmunoCAP FEIA method uses as the solid phase a flexible, hydrophobic cellulosic polymer to which allergen has been covalently linked. The advantage of this system is that it has a very high antigen binding capacity when compared to other systems and it has minimal non-specific binding with high total IgE.

**PDF Report**

No

**Day(s) and Time(s) Test Performed**

Monday through Friday

**Analytic Time**

1-2 days

**Maximum Laboratory Time**

3-6 days

**Performing Laboratory Location**

Eurofins Viracor

**Fees and Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

**CPT Code Information**

86003

**LOINC® Information**

Test ID	Test Order Name	Order LOINC Value
FDFEN	Dog Fennel ( <i>Anthemis cotula</i> ) IgE	7287-6

Result ID	Test Result Name	Result LOINC Value
Z5606	Dog Fennel ( <i>Anthemis cotula</i> ) IgE	7287-6
Z5607	Class	15705-7