
Overview**Method Name**

Spectrofluorometry (SF)

NY State Available

Yes

Specimen**Specimen Type**

Varies

Specimen Required

Submit only 1 of the following specimens:

Serum**Collection Container/Tube:** Red-top tube**Submission Container/Tube:** (T192) Amber plastic transport tube**Specimen Volume:** 1 mL**Collection Instructions:** Draw blood in a plain, red-top tube(s). **(Serum gel tube is not acceptable.)** Spin down and send 1 mL of serum frozen in an amber (T192) plastic transport tube to protect from light.**Note:** Label specimen appropriately (serum).**Plasma****Collection Container/Tube:** Lavender-top (EDTA)**Submission Container/Tube:** (T192) Amber plastic transport tube**Specimen Volume:** 1 mL**Collection Instructions:** Draw blood in a lavender-top (EDTA) tube(s). **(Plasma gel tube is not acceptable.)** Spin down and send 1 mL of EDTA plasma frozen in an amber (T192) plastic transport tube to protect from light.**Note:** Label specimen appropriately (plasma).**Specimen Minimum Volume**

0.3 mL

Reject Due To

Specimens other than	Serum, plasma
Hemolysis	NA
Lipemia	NA
Icteric	NA
Other	Not light protected; room temperature, refrigerated temperature, Polymer gel separation tubes (SST or PST)

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Frozen	7 days	LIGHT PROTECTED

Clinical and Interpretive

Reference Values

Reporting limit determined each analysis

Synonym(s): Dantrium

Usual therapeutic range: 0.2 - 3.5 mcg/mL

Performance

PDF Report

No

Day(s) and Time(s) Test Performed

Monday through Sunday

Analytic Time

7 day

Maximum Laboratory Time

9 - 11 days

Performing Laboratory Location

NMS Labs

Fees and Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80369

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FDANT	Dantrolene	9746-9

Result ID	Test Result Name	Result LOINC Value
Z1449	Dantrolene	9746-9
Z1859	Reporting Limit	19147-8