



# Test Definition: CALX

CALR Mutation Analysis, Myeloproliferative Neoplasm (MPN), Reflex, Varies

## Overview

### Useful For

Aiding in the distinction between a reactive cytosis and a chronic myeloproliferative disorder

Evaluating mutations in *CALR* in an algorithmic process for the MPNR / Myeloproliferative Neoplasm, *JAK2 V617F* with Reflex to *CALR* and *MPL*, Varies

### Method Name

Only orderable as a reflex. For more information see MPNR / Myeloproliferative Neoplasm, *JAK2 V617F* with Reflex to *CALR* and *MPL*, Varies.

Polymerase Chain Reaction (PCR) and Fragment Analysis

### NY State Available

Yes

## Specimen

### Specimen Type

Varies

### Specimen Required

Only orderable as a reflex. For more information see MPNR / Myeloproliferative Neoplasm, *JAK2 V617F* with Reflex to *CALR* and *MPL*, Varies.

**Submit only 1 of the following specimens:**

**Specimen Type:** Whole Blood

**Container/Tube:** Lavender top (EDTA) or yellow top (ACD solution B)

**Specimen Volume:** 3 mL

**Collection Instructions:**

1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. **Do not** aliquot.
3. Label specimen as blood.

**Specimen Stability Information:** Ambient (preferred)7 days/Refrigerate 7 days

**Specimen Type:** Bone marrow

**Container/Tube:** Lavender top (EDTA) or yellow top (ACD solution B)

**Specimen Volume:** 2 mL

**Collection Instructions:**

1. Invert several times to mix bone marrow.
2. Send bone marrow specimen in original tube. **Do not** aliquot.
3. Label specimen as bone marrow.

**Specimen Stability Information:** Ambient (preferred)7 days/Refrigerate 7 days

**Specimen Type:** Extracted DNA from blood or bone marrow

**Container/Tube:** 1.5- to 2-mL tube

**Specimen Volume:** Entire specimen

**Collection Instructions:**

1. Label specimen as extracted DNA from blood or bone marrow.
2. Provide volume and concentration of the DNA on the label.

**Specimen Stability Information:** Frozen (preferred)/Refrigerate/Ambient

**Specimen Minimum Volume**

Whole blood/Bone marrow: 1 mL

**Reject Due To**

Gross hemolysis	Reject
Paraffin shavings	Reject
Slides	Reject
Moderately to severely clotted	Reject
Paraffin embedded bone marrow aspirate clot or biopsy blocks	Reject

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Varies	7 days	

**Clinical & Interpretive**

**Clinical Information**

The *JAK2* (Janus kinase 2) gene codes for a tyrosine kinase (JAK2) associated with the cytoplasmic portion of a variety of transmembrane cytokine and growth factor receptors important for signal transduction in hematopoietic cells. Signaling

---

via JAK2 activation causes phosphorylation of downstream signal transducers and activators of transcription (STAT) proteins (eg, STAT5) ultimately leading to cell growth and differentiation. *BCR::ABL1*-negative myeloproliferative neoplasms (MPN) frequently harbor an acquired single nucleotide mutation in *JAK2* characterized as c.G1849T; p.Val617Phe (V617F). The *JAK2* V617F is present in 95% to 98% of polycythemia vera cases and 50% to 60% of primary myelofibrosis (PMF) and essential thrombocythemia (ET) cases. It has also been described infrequently in other myeloid neoplasms, including chronic myelomonocytic leukemia and myelodysplastic syndrome. Detection of the *JAK2* V617F is useful to help establish the diagnosis of MPN. However, a negative *JAK2* V617F result does not indicate the absence of MPN. Other important molecular markers in *BCR::ABL1*-negative MPN include *CALR* exon 9 mutation (20%-30% of PMF and ET) and *MPL* exon 10 mutation (5%-10% of PMF and 3%-5% of ET). Mutations in *JAK2*, *CALR*, and *MPL* are essentially mutually exclusive. A *CALR* mutation is associated with decreased risk of thrombosis in both ET and PMF and confers a favorable clinical outcome in PMF patients. A triple negative (*JAK2* V617F, *CALR*, and *MPL*-negative) genotype is considered a high-risk molecular signature in PMF.

### Reference Values

Only orderable as a reflex. For more information see MPNR / Myeloproliferative Neoplasm, *JAK2* V617F with Reflex to *CALR* and *MPL*, Varies.

An interpretive report will be provided.

### Interpretation

An interpretation will be provided under the MPNR / Myeloproliferative Neoplasm, *JAK2* V617F with Reflex to *CALR* and *MPL*, Varies.

### Cautions

A positive result is not specific for a particular subtype of myeloproliferative neoplasm and clinicopathologic correlation is necessary in all cases.

A negative result does not exclude the presence of a myeloproliferative neoplasm or other neoplastic process.

In rare cases, a mutation other than the V617F may be present in an area that interferes with primer or probe binding and cause a false-negative result.

### Clinical Reference

1. Klampfl T, Gisslinger H, Harutyunyan AS, et al. Somatic mutation of calreticulin in myeloproliferative neoplasms. *N Engl J Med*. 2013;369(25):2379-2390
2. Nangalia J, Massie CE, Baxter EJ, et al. Somatic *CALR* mutation in myeloproliferative neoplasms with nonmutated *JAK2*. *N Engl J Med*. 2013;369(25):2391-2405
3. Rumi E, Pietra D, Ferretti V, et al. *JAK2* or *CALR* mutation status defines subtypes of essential thrombocythemia with substantially different clinical course and outcomes. *Blood*. 2014;123(10):1544-1551
4. Rotunno G, Mannarelli C, Guglielmelli P, et al. Impact of calreticulin mutations on clinical and hematological phenotype and outcome in essential thrombocythemia. *Blood*. 2014;123(10):1552-1555
5. Tefferi A, Lasho TL, Finke CM, et al. *CALR* vs *JAK2* vs *MPL*-mutated or triple-negative myelofibrosis: clinical, cytogenetic and molecular comparisons. *Leukemia*. 2014;28(7):1472-1477
6. Greenfield G, McMullin MF, Mills K. Molecular pathogenesis of the myeloproliferative neoplasms. *J Hematol Oncol*. 2021;14(1):103

## Performance

### Method Description

Polymerase chain reaction (PCR) amplification of *CALR* exon 9 is performed on DNA isolated from the patient sample. The PCR product is then run on an ABI Genetic Analyzer for fragment analysis to detect insertions and deletions. An unmutated *CALR* will show an amplicon at 266 base pairs (bp), a mutated *CALR* with insertion will show an amplicon greater than 266 bp, and a mutated *CALR* with deletion will show an amplicon smaller than 266 bp. This assay has an analytical sensitivity of approximately 6% (ie, 6 mutation-containing cells in 100 total cells) in most mutation types, except for the rare type of 1-bp deletion, which has a sensitivity of approximately 20%. (Unpublished Mayo method)

### PDF Report

No

### Day(s) Performed

Monday through Friday

### Report Available

7 to 10 days

### Specimen Retention Time

Whole blood/Bone marrow: 2 weeks; Extracted DNA 3 months

### Performing Laboratory Location

Mayo Clinic Laboratories - Rochester Main Campus

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

### CPT Code Information

81219-CALR (calreticulin) (eg, myeloproliferative disorders), gene analysis, common variants in exon 9

### LOINC® Information

## Test Definition: CALX

CALR Mutation Analysis, Myeloproliferative  
Neoplasm (MPN), Reflex, Varies

Test ID	Test Order Name	Order LOINC® Value
CALX	CALR, Gene Mutation, Exon 9, Reflex	77174-1

Result ID	Test Result Name	Result LOINC® Value
36998	Final Diagnosis	22637-3