



Test Definition: HIVDR

HIV-1 Genotypic Drug Resistance to Reverse Transcriptase, Protease, and Integrase Inhibitors, Plasma

Overview

Useful For

Identifying HIV-1 genotypic mutations associated with resistance to nucleotide and non-nucleoside reverse-transcriptase inhibitors, protease inhibitors, and integrase strain transfer inhibitors

Guiding initiation or change of combination antiretroviral therapy in individuals, including children, with HIV-1 infection

Testing Algorithm

For information see [HIV Treatment Monitoring Algorithm](#)

Special Instructions

- [HIV Treatment Monitoring Algorithm](#)

Highlights

This assay uses next-generation sequencing to identify HIV-1 antiviral drug resistance-associated codon mutations in patients prior to or while receiving combination antiretroviral therapy. This test can be used to predict the likelihood of a favorable response to current US Food and Drug Administration-approved antiviral drug combinations used for treatment of HIV-1 infection.

Method Name

Reverse Transcription Polymerase Chain Reaction (RT-PCR) followed by Targeted Next-Generation Sequencing (NGS)

NY State Available

Yes

Specimen

Specimen Type

Plasma EDTA

Ordering Guidance

This test is intended for detection and identification of drug resistance-associated HIV-1 genotypic mutations in plasma specimens of individuals prior to or while receiving combination antiretroviral therapy.

Prior to requesting this test, patients must have a confirmed plasma HIV-1 RNA level (ie, viral load) of 1000 copies/mL or higher within the preceding 30 days. HIVQN / HIV-1 RNA Detection and Quantification, Plasma is available to provide this prerequisite test result. Alternately, if the patient's viral load is unknown, order HIQDR / HIV-1 RNA Quantification with Reflex to Genotypic Drug Resistance to Reverse Transcriptase, Protease, and Integrase Inhibitors, Plasma, which will

perform viral load followed by genotype, if appropriate.

For initial diagnosis of HIV, order HIVDX / HIV-1 and HIV-2 Antigen and Antibody Diagnostic Evaluation, Plasma.

Shipping Instructions

If shipment will be delayed for more than 24 hours, freeze plasma specimen at -70 degrees C (up to 60 days) until shipment on dry ice.

Necessary Information

The following ask-at-order entry question must be answered at the time of test ordering (mark answer on the test request form if not ordering electronically):

HIV-1 RNA level copies/mL in last 30 days = *(select answer option)*

- <1000
- 1000 to 1,000,000
- 1,000,001 to 10,000,000
- >10,000,000

Note: Test requests for submitted specimens with less than 1000 copies/mL (not sufficient amount for testing), "No," or no response entered will be canceled.

Specimen Required

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Lavender top (EDTA)

Submission Container/Tube: Plastic vial

Specimen Volume: 2.2 mL

Collection Instructions:

1. Centrifuge blood collection tube and aliquot plasma into plastic vial per collection tube manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).
2. Freeze aliquoted plasma for shipment.

Additional Information: Specimens submitted for HIV-1 genotyping must contain 1000 copies/mL or more of HIV-1 RNA.

Forms

If not ordering electronically, complete, print, and send a [Microbiology Test Request](#) (T244) with the specimen.

Specimen Minimum Volume

0.8 mL

Reject Due To

| | |
|-----------------|--------|
| Gross hemolysis | Reject |
| Gross lipemia | OK |

Test Definition: HIVDR

HIV-1 Genotypic Drug Resistance to Reverse Transcriptase, Protease, and Integrase Inhibitors, Plasma

| | |
|---------------|----|
| Gross icterus | OK |
|---------------|----|

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------|---------|-------------------|
| Plasma EDTA | Frozen (preferred) | 60 days | |
| | Refrigerated | 7 days | |

Clinical & Interpretive

Clinical Information

Antiviral resistance may compromise the efficacy of antiretroviral therapy (ART) in patients receiving such therapy for HIV1 infection. When combination therapy fails, detection and analysis of antiviral drug resistance-associated viral genotypic mutations can guide necessary changes to ART to suppress viral replication (ie, reduce viral load), thereby improving patient outcome.

HIV-1 is an RNA virus that infects cells and is then converted to complementary DNA by the action of the viral reverse transcriptase (RT) gene product. RT has little proofreading capacity and therefore, incorporates errors in the proviral DNA. These errors are transcribed into infectious viral particles when the proviral DNA is transcribed into RNA. Similarly, the enzyme protease catalyzes a polyprotein to produce peptides necessary for active viral replication. Although ART (combination of nucleoside and nonnucleoside reverse-transcriptase inhibitors, protease inhibitors, and/or integrase strain transfer inhibitors) may be effective in reducing the viral load, genotypic mutations arising in the drug-targeted HIV-1 genome due to selective pressure from antiviral therapy will result in antiviral resistance that may compromise such therapy.

Amplification and analysis of drug-targeted HIV-gene sequence allows identification of changes in nucleotide bases and associated amino acid codons that may cause antiviral drug resistance. Such genotypic changes are deemed as variants by comparing the sequence data of the patient's HIV strain to those of a wild-type HIV strain. The significance of these genotypic mutations in relation to antiviral resistance is then determined by a set of interpretive rules developed by a consensus panel of leading experts in the field of HIV-1 resistance. Relevant data presented at a recognized scientific conference or published in peer-reviewed journals are considered by the consensus panel in developing these rules. When necessary, reliable unpublished drug resistance data known to consensus panel members may be considered in the process. The interpretive rules are updated by the consensus panel annually after reviewing newly published data on HIV-1 genotypic drug resistance mutations.

Reference Values

An interpretive report will be provided.

Interpretation

Detectable HIV-1 genotypic mutations conferring resistance to an antiviral drug are reported as amino acid codon changes (eg, M184V) resulting from the nucleic acid base alterations, according to the interpretative algorithm of the Stanford HIV Database program. The codon mutations are categorized and interpreted in relation to previously

published data of phenotypic antiviral susceptibility tests on virus that harbor such mutations. Each codon mutation is assigned a drug penalty score. The total score generated from all mutations relevant to the specific antiviral drug is used to estimate the level of resistance to that drug. These interpretive rules may be updated periodically by the Stanford HIV Database Team after reviewing newly published data on HIV-1 genotypic drug resistance-associated codon mutations.

Susceptible (SUSC) indicates that the codon mutations present in patient's HIV-1 strain have not been associated with resistance to the specific drug (Stanford HIVdb total score 0 to 9).

Potential Low-Level Resistance (PLR) indicates that codon mutations detected have been associated with possible reduction in susceptibility to the specific drug (Stanford HIVdb score 10 to 14).

Low-Level Resistance (LR) indicates that codon mutations detected have been associated with reduction in susceptibility to the specific drug (Stanford HIVdb score 15 to 29).

Intermediate Resistance (IR) indicates that codon mutations detected have been associated with reduction in susceptibility to the specific drug (Stanford HIVdb score 30 to 59).

High-level Resistant (HR) indicates that codon mutations detected have been associated with maximum reduction in susceptibility to the specific drug (Stanford HIVdb > or = 60).

Unable to genotype indicates that the sequence data obtained are of poor quality to determine the presence or absence of resistance-associated codon mutations in the patient's HIV-1 strain. Probable causes of such poor sequence data include polymorphism in the region of the sequencing primers interfering with primer binding and subsequent sequencing reaction, or low viral load (ie, <1000 copies/mL).

Inconclusive indicates inability of the assay to reliably determine antiviral resistance because of the presence of polymerase chain reaction inhibitors or ambiguous or incomplete viral target sequences generated from the assay.

Cautions

Due to the complexity of the results generated, the International AIDS Society-USA Panel recommends expert interpretation of genotyping and phenotype test results for patient care management. A patient's response to antiviral therapy depends on multiple factors, including the percentage of patient's viral populations that is drug resistant, patient compliance with the prescribed drug therapy, patient access to adequate care, drug pharmacokinetics, and drug interactions. Drug resistance test results should be interpreted only in conjunction with clinical presentation and other laboratory markers when making therapeutic decisions.

Absence of resistance to a drug does not rule out the presence of reservoirs of drug-resistant virus in the infected individual.

The HIV-1 genotypic test is not a direct measure of drug resistance. Although genotypic testing can detect variants in the relevant HIV-1 genome, the significance of these variants requires careful interpretation to predict drug susceptibility. This assay's ability to amplify the target and detect genotypic mutation is poor and unreliable when the plasma HIV-1 viral load (VL) is less than 1000 copies/mL. Specimens submitted for this test should contain greater than or equal to

1000 copies/mL of HIV-1 RNA. Per assay manufacturer claims, the assay's ability to detect minor drug-resistant HIV-1 variants among 90% or more of HIV-1 group M strains varies depending on the VL in the tested plasma specimen; 20% or higher at VL of 1000 copies/mL, 10% or higher at VL of 5000 copies/mL, and 5% or higher at VL of 15,000 copies/mL.

The list of drug resistance-associated codon mutations and interpretive rules used by the Stanford HIV database are updated periodically by the Stanford HIV Database team. Therefore, the test results do not necessarily include all resistance-associated codon mutations described in the current medical literature.

Possible causes of treatment failure other than the development of drug resistance are poor adherence to medication regimen, drug potency, and individual variation in pharmacokinetics (eg, inadequate phosphorylation of nucleosides).

Clinical Reference

1. Weber J, Volkova I, Sahoo MK, Tzou PL, Shafer RW, Pinsky BA. Prospective evaluation of the Vela Diagnostics next-generation sequencing platform for HIV-1 genotypic resistance testing. *J Mol Diagn*. 2019;21(6): 961-970. doi:10.1016/j.jmoldx.2019.06.003
2. Avila-Rios S, Parkin N, Swanstrom R, et al. Next-generation sequencing for HIV drug resistance testing: laboratory, clinical, and implementation considerations. *Viruses*. 2020;12(6):617. doi: 10.3390/v12060617
3. Raymond S, Nicot F, Abravanel F, et al: Performance evaluation of the Vela Dx Sentosa next-generation sequencing system for HIV-1 DNA genotypic resistance. *J Clin Virol*. 2020;122:104229. doi:10.1016/j.jcv.2019.104229
4. Panel on Antiretroviral Guidelines for Adults and Adolescents: Guidelines for the use of antiretroviral agents in adults and adolescents with HIV. US Department of Health and Human Services. Updated September 12, 2024. Accessed March 27, 2025. Available at <https://clinicalinfo.hiv.gov/sites/default/files/guidelines/documents/adult-adolescent-arv/guidelines-adult-adolescent-arv.pdf>

Performance

Method Description

This test utilizes the US Food and Drug Administration-approved, commercially available Sentosa SQ HIV-1 Genotyping Assay, which is a next-generation sequencing assay based on a "sequencing by synthesis" method. The assay is designed to generate 2 amplicons (approximately 1500 base pairs [bp] and approximately 1000 bp in length) spanning the PR / RT- and INT-coding regions, respectively, of the HIV-1 genome for sequencing. Codon positions 1 to 99, 1 to 387, and 1 to 288 in the PR-, RT-, and INT-coding regions, respectively, are subsequently analyzed by the assay software for clinically relevant codon substitutions.

Clinical plasma specimens are subjected to automated HIV-1 RNA extraction and purification, followed by reverse-transcription polymerase chain reaction of HIV-1 target sequences, with both a system control and a positive control included in each assay run for quality control purposes. Automated DNA library preparation is performed using the amplified products, including enzymatic shearing, adapter ligation, purification, and normalization, prior to DNA template preparation and sequencing. Sequencing reactions are conducted with the Sentosa SQ301 sequencer, and the assembled sequence data are analyzed using proprietary analysis and interpretive software applications. HIV-1 antiviral

drug-resistance interpretations are based on algorithms implemented in the most current version of the Stanford University HIV Drug Resistance Database (HIVdb; Stanford University) using a 5% variant detection cutoff threshold set by the assay manufacturer. (Instruction manual: Sentosa SQ HIV-1 Genotyping Assay User Manual. Vela Diagnostics; Version 1.3, 05/2023)

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

3 to 10 days

Specimen Retention Time

60 days

Performing Laboratory Location

Mayo Clinic Laboratories - Rochester Superior Drive

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test has been cleared, approved, or is exempt by the US Food and Drug Administration and is used per manufacturer's instructions. Performance characteristics were verified by Mayo Clinic in a manner consistent with CLIA requirements.

CPT Code Information

0219U

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|------------------------------------|--------------------|
| HIVDR | HIV-1 Genotypic Drug Resistance, P | 90901-0 |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------------------------|---------------------|
| 616052 | HIV-1 Genotypic Drug Resistance, P | 80689-3 |

Test Definition: HIVDR

HIV-1 Genotypic Drug Resistance to Reverse Transcriptase, Protease, and Integrase Inhibitors, Plasma

| | | |
|--------|-------------------------------------|-----------------|
| 616729 | HIV-1 group M subtype | 100984-4 |
| 616737 | Nucleos(t)ide RT mutations | 45175-7 |
| 616918 | Reverse Transcriptase failed codons | 100983-6 |
| 616738 | Abacavir | 30287-7 |
| 616739 | Didanosine | 30284-4 |
| 616740 | Emtricitabine | 41402-9 |
| 616741 | Lamivudine | 30283-6 |
| 616742 | Stavudine | 30286-9 |
| 616743 | Tenofovir | 41396-3 |
| 616744 | Zidovudine | 30282-8 |
| 616745 | Nonnucleoside RT mutations | 45176-5 |
| 616746 | Doravirine | 91897-9 |
| 616747 | Efavirenz | 30291-9 |
| 616748 | Etravirine | 52749-9 |
| 616749 | Nevirapine | 30289-3 |
| 616750 | Rilpivirine | 68463-9 |
| 616751 | Protease Mutations | 33630-5 |
| 616919 | Protease failed codons | 100985-1 |
| 616752 | Atazanavir + Ritonavir | 49618-2 |
| 616753 | Darunavir + Ritonavir | 49630-7 |
| 616754 | Fosamprenavir + Ritonavir | 51409-1 |
| 616755 | Indinavir + Ritonavir | 49619-0 |
| 616756 | Lopinavir + Ritonavir | 42000-0 |
| 616757 | Nelfinavir | 30294-3 |
| 616758 | Saquinavir + Ritonavir | 49621-6 |
| 616759 | Tipranavir + Ritonavir | 49622-4 |
| 616760 | Integrase mutations | 61199-6 |
| 616920 | Integrase failed codons | 100986-9 |
| 616761 | Bictegravir | 90080-3 |
| 616762 | Cabotegravir | 96566-5 |
| 616763 | Dolutegravir | 72857-6 |
| 616764 | Elvitegravir | 72526-7 |
| 616765 | Raltegravir | 72525-9 |
| HIRVL | HIV RNA level copies/mL <30 days = | 89543-3 |
| 618206 | HIVDR_PR-RT_SEQ: | No LOINC Needed |
| 618207 | HIVDR_INT_SEQ: | No LOINC Needed |