

Overview

Method Name
High Performance Liquid Chromatography/ Tandem Mass Spectrometry (LC-MS/MS)

NY State Available
Yes

Specimen

Specimen Type
Serum Red

Specimen Required
Must submit one specimen per order. Specimens cannot be shared between multiple orders.
Specimen Type: Serum
Collection Container/Tube: Red top
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL
Collection Instructions:
1. Draw blood in a red-top tube(s). **Serum gel tube is not acceptable.**
2. Centrifuge and send 1 mL of serum refrigerated in a plastic, preservative-free vial.

Forms
If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

Specimen Minimum Volume
0.5 mL

Reject Due To

Other	Polymer gel separation tube (SST or PST)
-------	--

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	30 days	
	Ambient	14 days	
	Frozen		

Clinical & Interpretive

Clinical Information

Mean peak CBD plasma concentrations at 3-4 hours post-dose with Sativex at a low dose (5.4 mg of Delta-9 THC and 5.0 mg of Cannabidiol) were 1.6 +/- 0.4 ng/mL and at a high dose (16 mg of Delta-9 THC and 15 mg of Cannabidiol) were 6.7 +/- 2.0 ng/mL.

Following high dose 400 and 800 mg oral synthetic CBD in corn oil administration, mean peak CBD plasma concentrations occurred within 1.5-3 hours post-dose and were 181.2 +/- 39.8 and 221.1 +/- 35.6 ng/mL, respectively. The ratio of blood to serum or plasma concentration is unknown for this analyte.

Reference Values

Reporting limit determined each analysis

Units: ng/mL

Cautions

Known Interferences: Excessive amounts of oleamide from plastic tubes

Performance

PDF Report

No

Day(s) Performed

Monday-Sunday

Report Available

11 to 13 days

Performing Laboratory Location

NMS Labs

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80349
G0480, if appropriate

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FCBDS	Cannabidiol, Serum	82964-8

Result ID	Test Result Name	Result LOINC® Value
Z5547	Cannabidiol	82964-8