

Basic Metabolic Panel, Serum

Reporting Title: Basic Metabolic Panel, S

Performing Location: Rochester

Necessary Information:

Patient's age and sex are required.

Specimen Requirements:

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL Collection Instructions:

1. Serum gel tubes should be centrifuged within 2 hours of collection.

2. Red-top tubes should be centrifuged and the serum aliquoted into a plastic vial within 2 hours of collection.

Specimen Minimum Volume:

0.4 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated	24 hours	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
KS	Potassium, S	Numeric	mmol/L	2823-3
	Also used by tests: ELPSR, KS, NAKS			
NAS	Sodium, S	Numeric	mmol/L	2951-2
	Also used by tests: ELPSR, NAKS, NAS			
CL	Chloride, S	Numeric	mmol/L	2075-0
	Also used by tests: ELPSR, CL			
HCO3	Bicarbonate, S	Numeric	mmol/L	1963-8
	Also used by tests: HCO3, ELPSR			



Basic Metabolic Panel, Serum

Result ID	Reporting Name	Туре	Unit	LOINC®
AGAP	Anion Gap	Numeric		33037-3
	Also used by tests: ELPSR			
BUN	Bld Urea Nitrog (BUN), S	Numeric	mg/dL	3094-0
	Also used by tests: BUN			
CRTSA	Creatinine, S	Numeric	mg/dL	2160-0
	Also used by tests: CRTS1			
EGFR1	Estimated GFR (eGFR)	Numeric	mL/min/BSA	98979-8
	Also used by tests: CRTS1			
CA	Calcium, Total, S	Numeric	mg/dL	17861-6
	Also used by tests: CA			
GLURA	Glucose, Random, S	Numeric	mg/dL	2345-7
	Also used by tests: GLURA			

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
KS	Potassium, S			Yes	Yes
NAS	Sodium, S			Yes	Yes
CL	Chloride, S			Yes	Yes
HCO3	Bicarbonate, S			Yes	Yes
AGAP	Anion Gap			Yes	No
BUN	Bld Urea Nitrog (BUN), S			Yes	Yes
CRTS1	Creatinine with eGFR, S			Yes	Yes
CA	Calcium, Total, S			Yes	Yes
GLURA	Glucose, Random, S			Yes	Yes

CPT Code Information:

84132



Basic Metabolic Panel, Serum

84295

82435

82374

84520

82565

82310

82947

Reference Values:

SODIUM

<1 year: Not established > or =1 year: 135-145 mmol/L

POTASSIUM

<1 year: Not established > or =1 year: 3.6-5.2 mmol/L

CHLORIDE

<1 year: Not established 1-17 years: 102-112 mmol/L > or =18 years: 98-107 mmol/L

BICARBONATE

Males:

<1 year: Not established 1-2 years: 17-25 mmol/L 3 years: 18-26 mmol/L 4-5 years: 19-27 mmol/L 6-7 years: 20-28 mmol/L 8-17 years: 21-29 mmol/L > or =18 years: 22-29 mmol/L

Females:

<1 year: Not established 1-3 years: 18-25 mmol/L 4-5 years: 19-26 mmol/L 6-7 years: 20-27 mmol/L 8-9 years: 21-28 mmol/L > or =10 years: 22-29 mmol/L

ANION GAP

<7 years: Not established > or =7 years: 7-15

BLOOD UREA NITROGEN

Males:

<12 months: Not established 1-17 years: 7-20 mg/dL



Basic Metabolic Panel, Serum

> or =18 years: 8-24 mg/dL

Females:

<12 months: Not established 1-17 years: 7-20 mg/dL > or =18 years: 6-21 mg/dL

CREATININE

Males:

0-11 months: 0.17-0.42 mg/dL 1-5 years: 0.19-0.49 mg/dL 6-10 years: 0.26-0.61 mg/dL 11-14 years: 0.35-0.86 mg/dL > or =15 years: 0.74-1.35 mg/dL

Females:

0-11 months: 0.17-0.42 mg/dL 1-5 years: 0.19-0.49 mg/dL 6-10 years: 0.26-0.61 mg/dL 11-15 years: 0.35-0.86 mg/dL > or =16 years: 0.59-1.04 mg/dL

ESTIMATED GLOMERULAR FILTRATION RATE (eGFR)

>= 18 years old: > or =60 mL/min/BSA (body surface area)

Estimated GFR calculated using the 2021 CKD-EPI creatinine equation

Note: eGFR results will not be calculated for patients younger than 18 years old.

CALCIUM

<1 year: 8.7-11.0 mg/dL 1-17 years: 9.3-10.6 mg/dL 18-59 years: 8.6-10.0 mg/dL 60-90 years: 8.8-10.2 mg/dL >90 years: 8.2-9.6 mg/dL

GLUCOSE

0-11 months: Not established > or =1 year: 70-140 mg/dL