

Reporting Title: Renal Function Panel, S **Performing Location:** Rochester

Necessary Information:

Patient's age and sex are required.

Specimen Requirements:

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top Submission Container/Tube: Plastic vial Specimen Volume: 0.6 mL Collection Instructions: 1. Serum gel tubes should be centrifuged within 2 hours of collection. 2. Red-top tubes should be centrifuged and the serum aliquoted into a plastic vial within 2 hours of collection.

Specimen Minimum Volume:

0.5 mL

Forms:

If not ordering electronically, complete, print, and send a Renal Diagnostics Test Request (T830) with the specimen.

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------|----------|-------------------|
| Serum | Refrigerated | 24 hours | |

Result Codes:

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|--------------------------------------|---------|--------|--------|
| KS | Potassium, S | Numeric | mmol/L | 2823-3 |
| | Also used by tests: ELPSR, KS, NAKS | | | |
| NAS | Sodium, S | Numeric | mmol/L | 2951-2 |
| | Also used by tests: ELPSR, NAKS, NAS | | | |
| CL | Chloride, S | Numeric | mmol/L | 2075-0 |
| | Also used by tests: ELPSR, CL | | | |

Test Definition: RFAMA

Renal Function Panel, Serum

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|---------------------------------|---------|------------|---------|
| HCO3 | Bicarbonate, S | Numeric | mmol/L | 1963-8 |
| | Also used by tests: HCO3, ELPSR | | | |
| AGAP | Anion Gap | Numeric | | 33037-3 |
| | Also used by tests: ELPSR | | | |
| BUN | Bld Urea Nitrog (BUN), S | Numeric | mg/dL | 3094-0 |
| | Also used by tests: BUN | | | |
| CRTSA | Creatinine, S | Numeric | mg/dL | 2160-0 |
| | Also used by tests: CRTS1 | | | |
| EGFR1 | Estimated GFR (eGFR) | Numeric | mL/min/BSA | 98979-8 |
| | Also used by tests: CRTS1 | | | |
| CA | Calcium, Total, S | Numeric | mg/dL | 17861-6 |
| | Also used by tests: CA | | | |
| GLURA | Glucose, Random, S | Numeric | mg/dL | 2345-7 |
| | Also used by tests: GLURA | | | |
| ALB | Albumin, S | Numeric | g/dL | 1751-7 |
| | Also used by tests: LIVPR, ALB | | | |
| PHOS | Phosphorus (Inorganic), S | Numeric | mg/dL | 2777-1 |
| | Also used by tests: PHOS | | | |

LOINC and CPT codes are provided by the performing laboratory.

MAYO CLINIC LABORATORIES

Supplemental Report:

No

Components:

| Test ID | Reporting Name | CPT Units | CPT Code | Always Performed | Orderable Separately |
|---------|----------------|-----------|----------|---------------------|-------------------------|
| KS | Potassium, S | | | Yes | Yes |
| NAS | Sodium, S | | | Yes | Yes |
| CL | Chloride, S | | | Yes | Yes |
| HCO3 | Bicarbonate, S | | | Yes | Yes |

LABORATORIES

Test Definition: RFAMA

Renal Function Panel, Serum

| Test ID | Reporting Name | CPT Units | CPT Code | Always Performed | Orderable Separately |
|---------|---------------------------|-----------|----------|---------------------|-------------------------|
| AGAP | Anion Gap | | | Yes | No |
| BUN | Bld Urea Nitrog (BUN), S | | | Yes | Yes |
| CRTS1 | Creatinine with eGFR, S | | | Yes | Yes |
| CA | Calcium, Total, S | | | Yes | Yes |
| GLURA | Glucose, Random, S | | | Yes | Yes |
| ALB | Albumin, S | | | Yes | Yes |
| PHOS | Phosphorus (Inorganic), S | | | Yes | Yes |

CPT Code Information:

KS: 84132 NAS: 84295 CL: 82435 HCO3: 82374 BUN: 84520 CRTS1: 82565 CA: 82310 GLURA: 82947 ALB: 82040 PHOS: 84100

Reference Values:

SODIUM

<1 year: Not established > or =1 year: 135-145 mmol/L

POTASSIUM <1 year: Not established > or =1 year: 3.6-5.2 mmol/L

CHLORIDE <1 year: Not established 1-17 years: 102-112 mmol/L > or =18 years: 98-107 mmol/L

BICARBONATE Males <1 year: Not established 1-2 years: 17-25 mmol/L 3 years: 18-26 mmol/L 4-5 years: 19-27 mmol/L



Renal Function Panel, Serum

6-7 years: 20-28 mmol/L 8-17 years: 21-29 mmol/L > or =18 years: 22-29 mmol/L

Females

<1 year: Not established 1-3 years: 18-25 mmol/L 4-5 years: 19-26 mmol/L 6-7 years: 20-27 mmol/L 8-9 years: 21-28 mmol/L > or =10 years: 22-29 mmol/L

ANION GAP <7 years: Not established > or =7 years: 7-15

BLOOD UREA NITROGEN (BUN) Males <12 months: Not established 1-17 years: 7-20 mg/dL

Females <12 months: Not established 1-17 years: 7-20 mg/dL > or =18 years: 6-21 mg/dL

> or =18 years: 8-24 mg/dL

CREATININE

Males 0-11 months: 0.17-0.42 mg/dL 1-5 years: 0.19-0.49 mg/dL 6-10 years: 0.26-0.61 mg/dL 11-14 years: 0.35-0.86 mg/dL > or =15 years: 0.74-1.35 mg/dL

Females 0-11 months: 0.17-0.42 mg/dL 1-5 years: 0.19-0.49 mg/dL 6-10 years: 0.26-0.61 mg/dL 11-15 years: 0.35-0.86 mg/dL > or =16 years: 0.59-1.04 mg/dL

ESTIMATED GLOMERULAR FILTRATION RATE (eGFR)

>= 18 years old: > or =60 mL/min/BSA (body surface area) Estimated GFR calculated using the 2021 CKD-EPI creatinine equation. Note: eGFR results will not be calculated for patients younger than 18 years old.

CALCIUM <1 year: 8.7-11.0 mg/dL 1-17 years: 9.3-10.6 mg/dL



18-59 years: 8.6-10.0 mg/dL 60-90 years: 8.8-10.2 mg/dL >90 years: 8.2-9.6 mg/dL

GLUCOSE

0-11 months: Not established > or =1 year: 70-140 mg/dL

Total Protein > or =1 year: 6.3-7.9 g/dL Reference values have not been established for patients who are younger than 12 months of age.

ALBUMIN

> or =12 months: 3.5-5.0 g/dL
Reference values have not been established for patients who are younger than 12 months of age.

PHOSPHORUS

Males 1-4 years: 4.3-5.4 mg/dL 5-13 years: 3.7-5.4 mg/dL 14-15 years: 3.5-5.3 mg/dL 16-17 years: 3.1-4.7 mg/dL > or =18 years: 2.5-4.5 mg/dL Reference values have not been established for patients that are younger than 12 months of age.

Females 1-7 years: 4.3-5.4 mg/dL 8-13 years: 4.0-5.2 mg/dL 14-15 years: 3.5-4.9 mg/dL 16-17 years: 3.1-4.7 mg/dL > or =18 years: 2.5-4.5 mg/dL Reference values have not been established for patients that are younger than 12 months of age.