

### **Test Definition: CMAT**

Chromosomal Microarray, Tumor, Fresh or Frozen using Affymetrix Cytoscan HD

Reporting Title: Chromosomal Microarray, Tumor

Performing Location: Rochester

#### **Ordering Guidance:**

This test is not performed on formalin-fixed, paraffin-embedded (FFPE) specimens. If testing is needed for FFPE specimens, order CMAPT / Chromosomal Microarray, Tumor, Formalin-Fixed Paraffin-Embedded.

### **Necessary Information:**

- 1. A reason for testing must be provided for testing to be performed.
- 2. A pathology report should accompany the specimen. If this information is not available at the time of order, submit as soon as possible for appropriateness of testing and to aid in interpretation of results.

### **Specimen Requirements:**

Submit only 1 of the following specimens:

Supplies: Hank's Solution (T132)

Specimen Type: Tumor biopsy

Container/Tube: Sterile container with sterile Hank's balanced salt solution, Ringer's solution, or normal saline

Specimen Volume: 0.5-3 cm(3) or larger

Specimen Type: Lymph node

Container/Tube: Sterile container with sterile Hank's balanced salt solution, Ringer's solution, or normal saline.

Specimen Volume: 1 cm(3)

Specimen Type: Skin biopsy

Container/Tube: Sterile container with sterile Hank's balanced salt solution, Ringer's solution, or normal saline.

Specimen Volume: 4-mm diameter

Collection Instructions:

- 1. Wash biopsy site with an antiseptic soap.
- 2. Thoroughly rinse area with sterile water.
- 3. Do not use alcohol or iodine preparations.
- 4. A local anesthetic may be used.
- 5. Biopsy specimens are best taken by punch biopsy to include full thickness of dermis.

### Specimen Minimum Volume:

Tumor Biopsy: 3 cm(3) Lymph Node: 1 cm(3) Skin Biopsy: 4 mm diameter



## **Test Definition: CMAT**

Chromosomal Microarray, Tumor, Fresh or Frozen using Affymetrix Cytoscan HD

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Refrigerated		

## Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
CMAT	CG905	Reason for Referral	Plain Text	Yes

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
54729	Result	Alphanumeric		62356-1
54730	Nomenclature	Alphanumeric		62356-1
54728	Result Summary	Alphanumeric		50397-9
54731	Interpretation	Alphanumeric		69965-2
CG905	Reason for Referral	Alphanumeric		42349-1
54743	Specimen	Alphanumeric		31208-2
54732	Source	Alphanumeric		31208-2
54733	Method	Alphanumeric		85069-3
53424	Additional Information	Alphanumeric		48767-8
54734	Released By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

No

#### **CPT Code Information:**

81277



# **Test Definition: CMAT**

Chromosomal Microarray, Tumor, Fresh or Frozen using Affymetrix Cytoscan HD

### **Reference Values:**

An interpretive report will be provided.