

**Reporting Title:** Fish and Shellfish Panel IgG**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume:**

1 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	365 days	
	Ambient	7 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
Z3900	Clam IgG	Alphanumeric		60368-8
Z3901	Codfish/Scrod IgG	Alphanumeric		49099-5
Z3902	Crab IgG	Alphanumeric		60370-4
Z3903	Lobster IgG	Alphanumeric		60383-7
Z3904	Oyster IgG	Alphanumeric		60388-6
Z3905	Red Snapper IgG	Alphanumeric		60404-1
Z3906	Salmon IgG	Alphanumeric		60400-9
Z3907	Sardine/Pilchard IgG	Alphanumeric		63134-1
Z3908	Shrimp IgG	Alphanumeric		60403-3
Z3909	Sole IgG	Alphanumeric		60405-8
Z3910	Trout IgG	Alphanumeric		60409-0
Z3911	Tuna IgG	Alphanumeric		60410-8

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86001 x 12

**Reference Values:**

Reference ranges have not been established for food-specific IgG tests. The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for evaluation by diet elimination and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of the test.