

Reporting Title: HIV Ab Differentiation Prenatal, P**Performing Location:** Rochester**Ordering Guidance:**

This test is not suitable for follow-up testing of patients with reactive results from any rapid HIV tests, except for pregnant patients who were reactive by the Determine HIV-1/-2 Ag/Ab Combo rapid point-of-care test on serum or plasma (but not whole blood). Per the latest CDC recommended HIV testing algorithm patients with reactive results from any rapid HIV tests should be tested subsequently with laboratory-based HIV antigen and antibody combination immunoassays, such as HIVSP / HIV Antigen and Antibody Prenatal Routine Screen, Plasma.

If specimens are autopsy or cadaver blood sources, the proper US Food and Drug Administration licensed assay is HV1CD / HIV-1 and HIV-2 Antibodies for Cadaveric or Hemolyzed Specimens, Serum.

New York State clients: This test should not be requested for maternal/newborn HIV screening on specimens originating in New York State, due to state regulatory requirements for expedited result reporting.

Necessary Information:

Date of collection is required.

Specimen Requirements:

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Lavender top (EDTA)

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions:

1. Centrifuge blood collection tube per collection tube manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).
2. Aliquot plasma into plastic vial.

Specimen Minimum Volume:

0.8 mL

Specimen Type	Temperature	Time	Special Container
Plasma	Frozen (preferred)	30 days	
	Refrigerated	6 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
91952	HIV-1 Ab Differentiation Prenatal, P	Alphanumeric		68961-2
91953	HIV-2 Ab Differentiation Prenatal, P	Alphanumeric		81641-3

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86701
86702
87535 (if appropriate)
87538 (if appropriate)
87536 (if appropriate)

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
HPP12	HIV-1/HIV-2 RNA Detect Prenatal, P			No	Yes
HIQNP	HIV-1 RNA Detect/Quant Prenatal, P			No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
HPP12	616344	HIV-1 RNA	Alphanumeric		25835-0
HPP12	616345	HIV-2 RNA	Alphanumeric		69353-1
HIQNP	65567	HIV-1 RNA Detect/Quant Prenatal, P	Alphanumeric	copies/mL	70241-5

Reference Values:

Negative