

Thrombin Time (Bovine), Plasma

Reporting Title: Thrombin Time (Bovine), P **Performing Location:** Rochester

Specimen Requirements:

Specimen Type: Platelet-poor plasma Patient Preparation: Fasting preferred Collection Container/Tube: Light-blue top (3.2% sodium citrate) Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL Collection Instructions:

1. For complete instructions, see Coagulation Guidelines for Specimen Handling and Processing.

2. Centrifuge, transfer all plasma into a plastic vial, and centrifuge plasma again.

3. Aliquot plasma into a separate plastic vial leaving 0.25 mL in the bottom of centrifuged vial.

4. Freeze plasma immediately (no longer than 4 hours after collection) at -20 degrees C or, ideally, -40 degrees C or below.

Additional Information:

1. Double-centrifuged specimen is critical for accurate results as platelet contamination may cause spurious results.

2. Each coagulation assay requested should have its own vial.

Specimen Minimum Volume:

0.5 mL

Forms:

If not ordering electronically, complete, print, and send a Coagulation Test Request (T753) with the specimen.

| Specimen Type | Temperature | Time | Special Container | |
|---------------|-------------|---------|-------------------|--|
| Plasma Na Cit | Frozen | 14 days | | |

Result Codes:

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|---------------------------|---------|------|---------|
| TTSC | Thrombin Time (Bovine), P | Numeric | sec | 46717-5 |

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:



Thrombin Time (Bovine), Plasma

85670

Reference Values:

15.8-24.9 seconds