

Test Definition: MPLVS

MPL Exon 10 Mutation Detection, Varies

Reporting Title: MPL Exon 10 Mutation Detection, V **Performing Location:** Rochester

Specimen Requirements:

Submit only 1 of the following specimens:

Specimen Type: Peripheral blood
Container/Tube: Lavender top (EDTA) or yellow top (ACD)
Specimen Volume: 3 mL
Collections Instructions:
1. Invert several times to mix blood.
2. Send specimen in original tube. Do not aliquot.
3. Label specimen as blood.
Specimen Stability: Ambient (preferred)/Refrigerate

Specimen Type: Bone marrow Container/Tube: Lavender top (EDTA) or yellow top (ACD) Specimen Volume: 2 mL Collections Instructions: 1. Invert several times to mix bone marrow. 2. Send specimens in original tube. Do not aliquot.

3. Label specimen as bone marrow.

Specimen Stability: Ambient (preferred)/Refrigerate

Specimen Type: Extracted DNA from blood or bone marrow Container/Tube: 1.5- to 2- mL tube Specimen Volume: Entire specimen Collection Instructions: Label specimen as extracted DNA from blood or bone marrow and provide indication of volume and concentration of DNA. Specimen Stability Information: Frozen (preferred)/Refrigerated/Ambient

Specimen Minimum Volume:

Blood, Bone marrow: 0.5 mL; Extracted DNA: 50 mcL at 20 ng/mcL concentration

Forms:

If not ordering electronically, complete, print, and send a Hematopathology/Cytogenetics Test Request (T726) with the specimen.

| Specimen Type | Temperature | Time | Special Container |
|---------------|-------------|--------|-------------------|
| Varies | Varies | 7 days | |



MPL Exon 10 Mutation Detection, Varies

Ask at Order Entry (AOE) Questions:

| Test ID | Question ID | Description | Туре | Reportable |
|---------|-------------|---|-------------|------------|
| MPLVS | MP051 | Specimen Type: • Peripheral blood • Bone marrow • Extracted DNA from blood • Extracted DNA from bone marrow | Answer List | Yes |

Result Codes:

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|---------------------|--------------|------|---------|
| MP051 | Specimen Type | Alphanumeric | | 31208-2 |
| 602600 | Interpretation | Alphanumeric | | 69047-9 |
| 602601 | Signing Pathologist | Alphanumeric | | 19139-5 |

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

81339-MPL (myeloproliferative leukemia virus oncogene, thrombopoietin receptor, TPOR) (eg, myeloproliferative disorder), exon 10 sequence

Reference Values:

An interpretive report will be provided.