

## **Test Definition: NMRLP**

Nuclear Magnetic Resonance Lipoprotein Profile, Serum

Reporting Title: NMR Lipoprotein Profile, S

Performing Location: Rochester

#### **Specimen Requirements:**

Patient Preparation:

1. Fasting overnight (12-14 hours) is required. On night before examination, evening meal should be eaten before 6 p.m. and should contain no fatty foods.

2. Patient must not consume any alcohol for 24 hours before the specimen is collected.

Collection Container/Tube: Red top (serum gel/SST are not acceptable)

Submission Container/Tube: Plastic vial

Specimen Volume: 1.5 mL Collection Instructions:

1. Allow isopropyl alcohol (from phlebotomy site prep) to dry thoroughly before venipuncture.

2. Centrifuge and aliquot serum into a plastic vial.

#### **Specimen Minimum Volume:**

1 mL

#### Forms:

If not ordering electronically, complete, print, and send a Cardiovascular Test Request Form (T724) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Frozen	14 days	
	Ambient	8 hours	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
606167	LDL Particles, S	Alphanumeric	nmol/L	54434-6
606168	HDL Particles, S	Alphanumeric	mcmol/L	49748-7
606169	LDL Cholesterol (NMR), S	Alphanumeric	mg/dL	2089-1

LOINC and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

No



# **Test Definition: NMRLP**

Nuclear Magnetic Resonance Lipoprotein Profile, Serum

#### **CPT Code Information:**

83704

### **Reference Values:**

> or =18 years:

LDL Particles:

Desirable: <1,000 nmol/L

Above Desirable: 1,000-1,299 nmol/L Borderline high: 1,300-1,599 nmol/L

High: 1,600-2,000 nmol/L Very high: > or =2,000 nmol/L

HDL Particles: Male: >30 mcmol/L Female: >35 mcmol/L

LDL Cholesterol (NMR): Desirable: <100 mg/dL

Above Desirable: 100-129 mg/dL Borderline high: 130-159 mg/dL

High: 160-189 mg/dL Very high: > or =190 mg/dL

Reference values have not been established for patients younger than 18 years of age.