Hepatitis C Virus Genotypic Antiviral Drug Resistance, Serum

Test Definition: HCVDR

Reporting Title: HCV Genotypic Drug Resistance, S **Performing Location:** Rochester

Ordering Guidance:

This test is intended for detection of preexisting antiviral drug resistance-associated substitutions in individuals known to be infected with hepatitis C virus (HCV) genotype 1a, 1b, or 3 (any subtype) and being considered for HCV NS3, NS5A, and NS5B inhibitor combination therapy.

Additional Testing Requirements:

MAYO CLINIC

Prior to requesting this test, patients must have a confirmed serum or plasma hepatitis C virus (HCV) RNA level of 5000 IU/mL or higher within the preceding 30 days and a known HCV genotype result of 1a, 1b, or 3 (any subtype). The following tests are available to provide these prerequisite results:

-HCVQG / Hepatitis C Virus (HCV) RNA Quantification with Reflex to HCV Genotype, Serum

-HCVQN / Hepatitis C Virus (HCV) RNA Detection and Quantification by Real-Time Reverse Transcription-PCR, Serum -HCVG / Hepatitis C Virus Genotype, Serum

Shipping Instructions:

1. Ship specimen frozen on dry ice only.

2. If shipment will be delayed for more than 24 hours, freeze serum at -20 to -80 degrees C before shipment, and transport on dry ice.

Necessary Information:

The following 2 questions must be answered at the time of test ordering (if not ordering electronically, note the answers on the test request):

1. What Is the Hepatitis C Virus (HCV) RNA level in IU/mL within the last 30 days? Provide an answer using the following ranges:

- <5000

- 5000 to 1,000,000
- 1,000,001 to 10,000,000
- 10,000,001 to 100,000,000
- >100,000,000

Note: If the answer to this question is not answered or is "Unknown," testing will be canceled.

2. Does the patient have a known hepatitis C genotype of 1a, 1b, or 3 (any subtype)? Yes or No. Note: If the answer to this question is "No," testing will be canceled.

Specimen Requirements:

Supplies: Sarstedt Aliquot Tube, 5 mL (T914) Collection Container/Tube: Serum gel Submission Container/Tube: Plastic vial Specimen Volume: 2.5 mL



Collection Instructions:

1. Centrifuge blood collection tube per manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).

2. Aliquot serum into a plastic vial.

Specimen Minimum Volume:

1.6 mL

Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T732) with the specimen.

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------|---------|-------------------|
| Serum SST | Frozen (preferred) | 60 days | ALIQUOT TUBE |
| | Refrigerated | 7 days | ALIQUOT TUBE |

Ask at Order Entry (AOE) Questions:

| Test ID | Question ID | Description | Туре | Reportable |
|---------|-------------|---|-------------|------------|
| HCVDR | HCVLR | HCV RNA level IU/mL last 30 days =: • 5,000 to 1,000,000 • 1,000,001 to 10,000,000 • 10,000,001 to 100,000,000 • >100,000,000 | Answer List | Yes |
| HCVDR | HCVKG | Is known HCV genotype 1a, 1b, or 3?: • Yes • No | Answer List | Yes |

Result Codes:

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|----------------------------------|--------------|------|---------|
| 604410 | HCV Genotypic Drug Resistance, S | Alphanumeric | | 77202-0 |
| 604376 | HCV Genotype | Alphanumeric | | 92731-9 |
| 604377 | HCV NS3 RAS | Alphanumeric | | 73654-6 |
| 604378 | Sequence failure at codons: | Alphanumeric | | 92732-7 |
| 604379 | Glecaprevir resistance | Alphanumeric | | 92733-5 |

LABORATORIES

Test Definition: HCVDR

Hepatitis C Virus Genotypic Antiviral Drug Resistance, Serum

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|-------------------------------------|--------------|------|---------|
| 604380 | Grazoprevir resistance | Alphanumeric | | 82523-2 |
| 604382 | Voxilaprevir resistance | Alphanumeric | | 92734-3 |
| 604383 | HCV NS5A RAS | Alphanumeric | | 73655-3 |
| 604384 | Sequence failure at codons: | Alphanumeric | | 92732-7 |
| 604385 | Daclatasvir resistance | Alphanumeric | | 82379-9 |
| 604386 | Elbasvir resistance | Alphanumeric | | 82376-5 |
| 604387 | Ledipasvir resistance | Alphanumeric | | 82377-3 |
| 604388 | Pibrentasvir resistance | Alphanumeric | | 92735-0 |
| 604389 | Velpatasvir resistance | Alphanumeric | | 82520-8 |
| 604390 | HCV NS5B RAS | Alphanumeric | | 73655-3 |
| 604391 | Sequence failure at codons: | Alphanumeric | | 92732-7 |
| 604392 | Sofosbuvir resistance | Alphanumeric | | 82382-3 |
| 604393 | HCV Genotype | Alphanumeric | | 92731-9 |
| 604394 | HCV NS3 RAS | Alphanumeric | | 73654-6 |
| 604395 | Sequence failure at codons: | Alphanumeric | | 92732-7 |
| 604396 | Glecaprevir resistance | Alphanumeric | | 92733-5 |
| 604397 | Grazoprevir resistance | Alphanumeric | | 82523-2 |
| 604399 | Voxilaprevir resistance | Alphanumeric | | 92734-3 |
| 604400 | HCV NS5A RAS | Alphanumeric | | 73655-3 |
| 604401 | Sequence failure at codons: | Alphanumeric | | 92732-7 |
| 604402 | Daclatasvir resistance | Alphanumeric | | 82379-9 |
| 604403 | Elbasvir resistance | Alphanumeric | | 82376-5 |
| 604404 | Ledipasvir resistance | Alphanumeric | | 82377-3 |
| 604405 | Pibrentasvir resistance | Alphanumeric | | 92735-0 |
| 604406 | Velpatasvir resistance | Alphanumeric | | 82520-8 |
| 604407 | HCV NS5B RAS | Alphanumeric | | 73655-3 |
| 604408 | Sequence failure at codons: | Alphanumeric | | 92732-7 |
| 604409 | Sofosbuvir resistance | Alphanumeric | | 82382-3 |
| HCVLR | HCV RNA level IU/mL last 30 days = | Alphanumeric | | 86955-2 |
| HCVKG | Is known HCV genotype 1a, 1b, or 3? | Alphanumeric | | 86955-2 |

LOINC and CPT codes are provided by the performing laboratory.



Supplemental Report:

No

CPT Code Information:

87900 87902 87999 (if appropriate for government payers)

Reference Values:

An interpretive report will be provided.